




+974 50066887
info@dawamedical.qa
45, Street 445, Zone 71
Doha, Qatar



Disclaimer

All product names, logos, and brands are the property of their respective owners. Dawa Medical showcases these items for informational and distribution purposes only, with full respect to the original manufacturers. 

جميع أسماء المنتجات والشعارات والعلامات التجارية هي ملك لأصحابها الأصليين. تقوم شركة دواء ميديكال بعرض هذه المنتجات لأغراض إعلامية وتوزيعية فقط، مع كامل الاحترام للشركات المصنعة الأصلية.



ABOUT BENECARE Since our establishment in 2001, BeneCare™ Medical (a trading name of Benefoot UK Ltd) has been a leading supplier of immobilisation and rehabilitation products to the National Health Service and Private Sector Hospitals. Our ranges include Casting Materials, Low Temperature Thermoplastics, Upper + Lower Limb Orthotics and Orthopaedic Footwear.



WHY US? We are a family owned company with a highly experienced team, dedicated to offering unbeatable service and long term relationships with all of our customers both in the UK and Worldwide. As a company we pride ourselves not only on our high level of personal service, but by sourcing our products from some of the finest manufacturers around the world. We can ensure that our quality is always of the highest standard whilst keeping our prices competitive and keen.

LOCATION + EXPANSION Our centrally situated warehouse and offices in the heart of Manchester ensures that we have daily pickups and collections. All orders received before 2pm are dispatched the same day on overnight carriers.

With regards to expansion, we are always looking to expand further and continue to provide you with products that improve patient wellbeing and help reduce clinic times. If you can't find what you need within our catalogue, please feel free to contact the office directly or your personal sales executive.



OUR MISSION BeneCare's overall mission is to develop and provide innovative and superior specialist products, services and education, which help to accelerate the rehabilitation of users and are universally accessible to all.



CONTENTS

CASTING & SPLINTING

04 - 09

ORTHOPAEDIC FOOTWEAR

10 - 16

ANKLE BRACES

18 - 19

KNEE SUPPORTS

20 - 23

ARM & SHOULDER

24 - 26

WRIST & THUMB

28 - 31

BENEKIDZ

32 - 36

PRODUCT INFO & INDEX

38 - 46

ORDERING INFORMATION

47

OUR VISION Our vision is to be the preferred medical supplies specialist, leading the industry in terms of quality, productivity and unrivalled customer care and support that help to maintain and restore personal independence.



info@dawamedical.qa
+974 50066887

CASTING MATERIALS

BeneCare have developed a range of casting materials to meet the needs of today's diverse treatment regimes. Developed in conjunction with leading healthcare professionals to ensure the exacting standards are met for today's NHS.

The range includes 3 polyester materials and one fibreglass, the versatility of the material will enable the clinician to choose the most suitable material for the indication.

BeneCast Advanced - Polyester casting tape benefits from a four way stretch allowing superior conformability, excellent lamination, fast setting, a smooth outer surface and soft edges. Suitable for FRC.

BeneCast Fibre Tx – Fibreglass casting tape with excellent lamination, great conformability, fast setting and low profile.

BeneCast Select – Polyester casting tape, the benefits include fast setting for a polyester, good lamination, excellent durability, and smooth surface. Benecast Flex – a Semi Rigid Polyester casting tape, excellent lamination, smooth surface, soft edges, suitable to create a cast of varying levels of rigidity. Removal with scissors or unwinding.

The versatility of these casting materials are the preferred choice for a variety of healthcare professionals, these include Orthopaedic Practitioners, Podiatrists, Paediatric and Neuro Physiotherapists.

The techniques these materials are suitable for include Conventional Casting, Focussed Rigidity Casting (FRC), Flexible Stabilisation (Combi casting), Functional Fracture Bracing. Paediatrics, Podiatry Off Loading DFU, Neuro splinting. The BeneCast casting materials are lightweight, breathable and provide ease of use for the end user, they help improve quality care for the patient, with promoting rehabilitation to a maximum attainable degree.

ADVANCED CASTING TAPE

05

FIBREGLASS CASTING TAPE

05

SELECT CASTING TAPE

06

FLEX CASTING TAPE

06

PLASTER OF PARIS

07

UNDERCAST PADDING

07

SPLINT ON A ROLL

07



PRE CUT MOULDED SPLINTS

08

COTTON CREPE BANDAGE

08

COTTON RIBBED STOCKINETTE

08

NON-WOVEN STERILE SWABS

09

NON-WOVEN SWABS

09

PRE-TAPE

09

BENECARE ADVANCED POLYESTER CASTING TAPE

BeneCast™ Advanced is our premium polyester casting material, made with a unique 4 way stretch, it benefits from a higher level of control and conformability, which provides the clinician with exacting characteristics to provide excellent comfort and support for the patient in modern fracture management and soft tissue indications.

BeneCast™ Advanced's key functions are to provide a casting material suitable for conventional casting, in addition to having clinical evidence to support the use of FRC (Focussed Rigidity Casting) using varying levels of rigidity appropriate to the indication.

The low profile breathable material benefits from soft edges, a smooth surface and fits to the contours of the limb perfectly, providing optimum stability to the limb.



FEATURES & BENEFITS

- Unique 4 Way stretch.
- Excellent handling characteristics.
- Excellent lamination with low profile and smooth surface.
- Clinical evidence to suggest suitable for FRC.
- Superior conformability.
- Lightweight and strong – provides maximum strength to weight ratio.
- Air permeability – allows evaporation of moisture from the skin.
- Excellent X-Ray transparency.

INDICATIONS OF USE

- Orthopaedic and trauma conditions.
- Fracture Management.
- Soft tissue injuries.
- DFU.
- Neurological symptoms.
- Paediatric indications.

BENECARE CODE: CCT

PLEASE TURN TO PAGE 38 FOR SIZING & FITTING OPTIONS

BENECARE FIBRE TX CASTING TAPE



BeneCast™ Fibre TX, is a synthetic fibre based casting tape that provides high quality patient care and comfort. By combining a fibreglass substrate with an advanced water activated polyurethane resin system.

BeneCast™ Fibre TX boasts a low tack application with a smooth finish and excellent lamination, providing a strong durable cast. The special air permeability properties also allows moisture to evaporate from the skin helping to reduce irritation.

FEATURES & BENEFITS

- Lightweight and strong - provides excellent strength to weight ratio.
- Air permeability - allows evaporation of moisture from the skin.
- Short application time.
- Smooth finish.
- Low profile.
- Fast setting.

INDICATIONS OF USE

- Orthopaedic and Trauma Indications.
- Fracture Management.
- Soft Tissue injuries.
- DFU.
- Paediatric indications.

BENECARE CODE: FCT

PLEASE TURN TO PAGE 38 FOR SIZING & FITTING OPTIONS

BENECAST SELECT CASTING TAPE



BeneCast™ Select, polyester casting tape, has excellent handling characteristics, low profile with a smooth surface, with exceptional durability. Ideal for primary and secondary casting, providing strength and varying levels of rigidity for maximum stability of the injury.

It also provides a lightweight and breathable cast which is easy to apply and conforms to body contours for easier cast application.

FEATURES & BENEFITS

- Lightweight but strong - provides excellent strength to weight ratio.
- Breathable therefore prevents maceration of the skin.
- 100% X-Ray transparency.
- Short application time.
- Low profile.
- Fast Setting.
- Excellent durability.

INDICATIONS OF USE

- Orthopaedic and trauma conditions.
- Fracture Management.
- Soft tissue injuries.
- DFU.
- Neurological symptoms.
- Paediatric indications.

BENECARE CODE:

CT

PLEASE TURN TO PAGE 38 FOR SIZING & FITTING OPTIONS

BENECAST FLEX CASTING TAPE

BeneCast™ Flex, semi-rigid casting tape, is polyester based with a polyurethane resin. The composition of the material provides semi-rigid support for the injury when used in isolation, and maximum stability when used in association with a rigid casting material to produce a combicast, whilst allowing functionality to the soft tissues.

It has excellent handling characteristics with ease of application for the end user conforming to all body contours.

Flex can be applied wet or dry, and removed with scissors or by unwinding.



FEATURES & BENEFITS

- BeneFlex™ is easy to use.
- BeneFlex™ is tack-free and clean to apply.
- High porosity and good moisture resistance allows BeneFlex™ casts to breath and dry out easily.
- Reduces risk of body odour and skin maceration under the cast.
- Easy to remove either by unwrapping or cutting with bandage scissors.
- Apply wet or dry.
- Ideal for paediatric patients.
- Flexible stabilisation.
- Low profile.
- Quick setting.
- Ease of application.

INDICATIONS OF USE

- Orthopaedic and trauma conditions.
- Modern closed fracture management.
- Soft Tissue injuries.
- Neuro splinting.
- Paediatric indications.
- DFU.
- Therapeutic casting and splinting.

BENECARE CODE:

CTSR

PLEASE TURN TO PAGE 39 FOR SIZING & FITTING OPTIONS

BENECARE PLASTER OF PARIS



BeneCast™ POP combines all the technician/patient benefits of the finest leno-weave gauze. Impregnated with a specially formulated blend of alpha and beta Plaster of Paris crystals incorporating high-grade binders and accelerants. It is then spooled onto a rigid plastic core or layered into slabs.

FEATURES & BENEFITS

- Excellent wetting-out of the bandage upon immersion into water.
- Excellent moulding to the body contours.
- Fast setting, strong and durable.
- Suitable for all trauma and orthopaedic conditions.

INDICATIONS OF USE

- Orthopaedic and Trauma conditions.
- Primary care casting.
- A&E.
- Fracture management.

BENECARE CODE: P.O.P

PLEASE TURN TO PAGE 39 FOR SIZING & FITTING OPTIONS

BENECARE UNDERCAST PADDING

Benecare Undercast Padding is a synthetic and protective material which blends excellent conformability and good cushioning with high patient care. The padding can be torn easily and has excellent conformability to the contours of the limbs, minimising patient trauma and enhancing ease of application.

FEATURES & BENEFITS

- Highly conformable, shapes easily and readily to all body contours and joints.
- Easily tearable, it has optimum strength to ensure that it will tear apart with no pain or discomfort to the patient.
- Feathers out neatly, clings to itself and blends in to prevent loops or ridges, to create a smooth finish.
- Low moisture absorbency, minimising skin maceration.
- Quick and easy to apply.

INDICATIONS OF USE

- For comfort and protection under a cast.
- Added compression under a bandage.
- For protection of a wound.



BENECARE CODE: SUP

PLEASE TURN TO PAGE 39 FOR SIZING & FITTING OPTIONS

BENECARE SPLINT ON A ROLL



BeneCast™ Splint on a Roll, Pre padded fibreglass splinting system benefits from a water repellent polypropylene felt padding. It is packed in a re-sealable foil pouch which minimizes waste, enhances productivity and gives the clinician the ability to cut the splint to the required length with no mess.

This versatile material has excellent conformability and a great alternative to Plaster of Paris. It is lightweight, strong and durable and can be removed easily by the patient or clinician for personal hygiene, woundcare and physiotherapy where clinically acceptable.

FEATURES & BENEFITS

- Fast and efficient splinting system.
- Strong & durable.
- Excellent quality polypropylene water repellent felt padding.
- Excellent conformability.
- Easy to use resealable foil pouch.
- Setting time 3-7 minutes with functional strength in 30 minutes.
- Splint can be washed in warm water with a thorough rinsing.

INDICATIONS OF USE

- Trauma & fracture management.
- Soft Tissue splinting.
- Neuro resting splint.

BENECARE CODE: SR

PLEASE TURN TO PAGE 39 FOR SIZING & FITTING OPTIONS

BENECARE PRE CUT MOULDED SPLINTS

Strong synthetic splints with water repellent felt padding. Designed for ease of use with a variety of splints available.

For upper and lower limb. Ideal for A & E, Theatres and Ambulance Services.



FEATURES & BENEFITS

- Can be used quickly in emergency situations by ambulance services.
- 100% penetration of X-ray.
- 5 x lighter & stronger than Plaster Casting.
- Fast/efficient way to make a removable splint.
- Secure using Benecare cotton crepe bandage.
- Pre-cut to its individual shape.
- Requires minimal trimming or preparation.
- Excellent conformability for shaping to the contours of the limb.

INDICATIONS OF USE

- Trauma & orthopaedics.
- In situations where ease of use is crucial.

BENECARE CODE:

PCMS

PLEASE TURN TO PAGE 39 FOR SIZING & FITTING OPTIONS

BENECARE COTTON CREPE BANDAGE

BeneCare Cotton Crepe Bandage conforms easily to all body contours and is a first aid essential. It can be used to secure dressings in place over wounds, to support the limb following a strain, sprain or simple fracture. Will provide excellent compression to help control the occurrence of oedema.



FEATURES & BENEFITS

- Provides padding/protection and support.
- Exceptionally high stretch.
- Excellent conformability to the contours of the limb.
- Soft and supple.
- Exceptional quality, available in a variety of sizes.
- Permeable to air and unaffected by ointments and perspiration.
- Extremely durable high twist cotton construction.

INDICATIONS OF USE

- To provide compression for sprains and strains.
- Secures splints and POP back slabs.
- Secures dressings.

BENECARE CODE:

SCCB

PLEASE TURN TO PAGE 40 FOR SIZING & FITTING OPTIONS

COTTON RIBBED STOCKINETTE

BeneCare Cotton Ribbed Stockinette is an inexpensive and convenient form of covering which makes it especially suitable for use as an undercast stockinette.



FEATURES & BENEFITS

- Protects skin from friction under synthetic and plaster casts.
- 100% cotton fabric is comfortable on the patients skin.
- Knitted construction allows air to circulate, helping avoid maceration of protected areas.
- The conformability provides protection and comfort for the patient.

INDICATIONS OF USE

- Under a cast.
- To secure dressings.
- Protects bedding and clothing from ointments and creams.

BENECARE CODE:

SN

PLEASE TURN TO PAGE 40 FOR SIZING & FITTING OPTIONS

NON WOVEN STERILE SWABS



BeneCare™ Non Woven Sterile Swabs 7.5 x 7.5cm packed in peel-apart pouches 5 swabs per pouch, 20 per box. The non woven sterile swabs are extremely soft and gentle to the touch yet strong and durable when wet. Manufactured to comply with the Drug Tariff Technical Specification 28, these swabs have high levels of absorption, fast wicking and are extremely low linting.

BENECARE CODE: NWSS

PLEASE TURN TO PAGE 40 FOR SIZING & FITTING OPTIONS

FEATURES & BENEFITS

- Low adhesion, important when treating open wounds.
- Sterile - in peel apart pouches in 5 pieces.
- 20 pouches per box.
- Conforms to Drug Tariff Specification 28.

INDICATIONS OF USE

- Suitable for an indication where high levels of absorption are required.

NON WOVEN SWABS

The 4 ply non woven swab is made from top grade absorbent fabrics. They have a standard 70% Viscose, 30% Polyester fibre mix, with a 30gsm weight. They provide improved strength, absorbency and scrubbing texture compared to traditional cotton gauze swabs. Manufactured to comply with the Drug Tariff Technical Specification 28



FEATURES & BENEFITS

- Higher absorbency than 12-ply cotton gauze swabs.
- Made from 4 ply woven top grade absorbent fabric.
- Less linting than traditional gauze.
- Holds its shape and texture better than cotton gauze swabs.
- Reduces snagging on structures and provides excellent wound release.
- Versatile use, cleans, absorbs and wicks.
- Conforms to Drug Tariff Specification 28.

INDICATIONS OF USE

- Suitable for an indication where high levels of absorption are required.

BENECARE CODE: NWS

PLEASE TURN TO PAGE 40 FOR SIZING & FITTING OPTIONS

BENECARE PRE TAPE



BeneCare Pre-Tape is a thin tape for applying over padding to minimise ridging and create a smooth compressed layer underneath Plaster of Paris and synthetic casts. BeneCare Pre-Tape underwrap is a porous and water resistant material which allows the skin to breathe, which helps prevent maceration.

FEATURES & BENEFITS

- Particularly useful when creating windows for inspection and for easy conversion of casts to resting splints.
- Latex free.
- Added patient comfort.
- Improves removal of casts.

INDICATIONS OF USE

- Under a cast to prevent ridging.
- To enable a bivalve or window to be removed with ease.

BENECARE CODE: BPT

PLEASE TURN TO PAGE 40 FOR SIZING & FITTING OPTIONS

ORTHOPAEDIC FOOTWEAR



The need for a medical/orthotic shoe is ever increasing, with the diversity of lifestyle choices and activities, this would include both older and younger generations.

BeneCare have a comprehensive high quality range of orthopaedic shoes which are suitable for a variety of indications and requirements, ranging from the Original Post Op and High Top Shoe, Square Toe Medical shoe and the Wedge shoe.

This diverse range will be the ideal choice for a number of indications, including general foot situations where a regular shoe does not fit over a cast, diabetic foot ulcers, post-surgical procedures of the foot, and chronic inflammatory disorders.

ORIGINAL POST OP MEDICAL SHOE 11

POST-OP SHOE WITH TOE COVER 11

HIGH TOP MEDICAL SHOE 12

HIGH-TOP SHOE WITH TOE COVER 12

SQUARE TOE SHOE 12

SQUARE TOE SHOE WITH TOE COVER 13

WEDGE SHOE 13

PAEDIATRIC CAST SHOE 13

B-COMFY SHOE 14



B-COMFY BOOT 14

ECO WALKER 14

FIXED WALKER 15

ROM WALKER 15

BENEKIDZ WALKER 15

SWEDISH AFO 16

FPFS SUPPORT 16

ACHILLES HEEL PADS 16



BENEFOOT ORIGINAL POST-OP MEDICAL SHOE

The BeneFoot Original Post-Op Medical shoe benefits from a rounded toe design, it has a strong rigid sole with a slight rocker bottom providing excellent support and assists with ambulation. The deep tread will prevent slipping by creating good friction on hard surfaces. The black nylon mesh body of the shoe has a padded collar with a reinforced heel counter for added comfort and protection and the dual density insole provides comfort and allows modification when required. The extra long straps provide an adjustable secure hook and loop closure system that can be adjusted easily over dressings. The rigid sole of the Original Post Op shoe has the capabilities to distribute the plantar pressure and also help prevent any secondary musculoskeletal injury from occurring. This shoe is therefore the ideal choice for the treatment of diabetic foot ulcers where off-loading is required.



FEATURES & BENEFITS

- Strong rigid sole with excellent tread.
- Slight rocker bottom.
- Breathable nylon mesh body.
- Padded collar and reinforced heel counter.
- Dual density insole.
- Adjustable hook and loop closure system.


INDICATIONS OF USE

- Post Operative.
- Over large dressings and bandages.
- Diabetic foot conditions including plantar wounds.
- Over a lower limb cast.
- Any situation where the size of the foot prohibits normal footwear.

BENECARE CODE:

PO

PLEASE TURN TO PAGE 40 FOR SIZING & FITTING OPTIONS



BENEFOOT ORIGINAL POST-OP MEDICAL SHOE (WITH TOE COVER)

Benefits from the rounded toe design with toe cover. A strong rigid sole with a slight rocker bottom, provides excellent support to assist with ambulation.

The deep tread will prevent slipping by creating good friction on hard surfaces. The black nylon mesh body of the shoe has a padded collar with a reinforced heel counter for added comfort and protection.

The dual density insole provides excellent comfort and allows modification when required. The extra long straps provide a secure hook and loop closure system that can be adjusted easily over dressings etc.



FEATURES & BENEFITS

- Strong rigid sole with excellent tread.
- Slight rocker bottom.
- Breathable nylon mesh body.
- Padded collar and reinforced heel counter.
- Dual density insole.
- Adjustable hook and loop closure system.

INDICATIONS OF USE

- Post Operative.
- Over large dressings and bandages.
- Diabetic foot conditions including plantar wounds.
- Over a lower limb cast.
- Any situation where the size of the foot prohibits normal footwear.

BENECARE CODE:

POTC

PLEASE TURN TO PAGE 40 FOR SIZING & FITTING OPTIONS





BENECARE CODE: POHT

PLEASE TURN TO PAGE 40 FOR SIZING & FITTING OPTIONS

BENEFoot HIGH TOP MEDICAL SHOE

The BeneFoot High Top Medical Shoe benefits from a rounded toe design. It has a strong rigid sole with a slight rocker bottom providing excellent support to assist with ambulation. The deep tread will prevent slipping by creating good friction on hard surfaces. The black nylon mesh body of the shoe has added ankle support with a reinforced heel counter for added comfort and protection.

The extra ankle strap will prevent slippage and forward foot movement, consequently reducing shearing forces of the foot and the dual density insole provides more comfort and allows modification when required. The extra long straps provide an adjustable secure hook and loop closure system that can be adjusted easily over dressings, oedematous legs or a cast.

FEATURES & BENEFITS

- Strong rigid sole with excellent tread.
- Slight rocker bottom.
- Breathable nylon mesh body.
- Ankle strap eliminates heel slippage and forward foot movement.
- Extended ankle support & reinforced heel counter.
- Dual Density insole.
- Adjustable hook and loop closure system.

INDICATIONS OF USE

- Post operative.
- Over large dressings and bandages.
- Diabetic foot conditions including plantar wounds.
- Over a lower limb cast.
- Any situation where the size of the foot prohibits normal footwear.

BENEFoot HIGH TOP MEDICAL SHOE (WITH TOE COVER)

Like the 'Original' this shoe benefits from a rounded toe with toe cover, a strong rigid sole with the slight rocker bottom. It provides excellent support that assists with ambulation, and the deep tread prevents slipping. The black nylon mesh body of the shoe has added ankle support with a reinforced heel counter for added comfort and protection. The extra ankle strap will prevent slippage and forward foot movement, consequently reducing shearing forces of the foot. The dual density insole provides more comfort and allows modification when required.



BENECARE CODE: POHT TC

PLEASE TURN TO PAGE 40 FOR SIZING & FITTING OPTIONS

FEATURES & BENEFITS

- Strong rigid sole with excellent tread.
- Slight rocker bottom.
- Breathable nylon mesh body.
- Ankle strap eliminates heel slippage and forward foot movement.
- Extended ankle support and reinforced heel counter.
- Dual Density insole.
- Adjustable hook and loop closure system.

INDICATIONS OF USE

- Post operative.
- Over large dressings and bandages.
- Diabetic foot conditions including plantar wounds.
- Over a lower limb cast.
- Any situation where the size of the foot prohibits normal footwear.



BENECARE CODE: BMST

PLEASE TURN TO PAGE 40 FOR SIZING & FITTING OPTIONS

BENEFoot MEDICAL SQUARE TOE SHOE

The BeneFoot Medical Square Toe Shoe has a square toe design that provides additional protection for the wider foot and the toes when K-wires are present. It benefits from adjustable strapless closure system on the forefoot which helps eliminate buckle pressure, with an additional ankle strap to prevent forefoot slippage.

The strong rigid sole, with a slight rocker bottom, assists with a more natural gait whilst mobilising.

FEATURES & BENEFITS

- Square toe design.
- Slight rocker bottom sole aids ambulation.
- Ankle strap eliminates heel slippage.
- Strapless forefoot closure system.
- Removable insole fully customisable.

INDICATIONS OF USE

- Post operative.
- Diabetic foot conditions, including the peg assist insole.
- To accommodate large dressings.
- Any situation where the size of foot prohibits normal footwear.
- Can be used over a cast.

BENEFoot MEDICAL SQUARE SHOE (WITH TOE COVER)

The BeneFoot Medical Square Toe Shoe (With Toe Cover) has a square toe design that provides additional protection for the wider foot and the toes when K-wires are present. It benefits from adjustable strapless closure system on the forefoot which helps eliminate buckle pressure, with an additional ankle strap to prevent forefoot slippage.

The strong rigid sole with a slight rocker bottom assists with a more natural gait whilst mobilising.

FEATURES & BENEFITS

- Square toe design.
- Slight Rocker bottom sole aids ambulation.
- Ankle strap eliminates heel slippage.
- Strapless forefoot closure system.
- Removable insole fully customisable

INDICATIONS OF USE

- Post operative.
- Diabetic foot conditions including the peg assist insole.
- To accommodate large dressings.
- Any situation where the size of foot prohibits normal footwear.
- Can be used over a cast.



BENECARE CODE:

BMST TC

PLEASE TURN TO PAGE 41 FOR SIZING & FITTING OPTIONS

BENEFoot WEDGE SHOE

The BeneFoot Wedge Shoe has been clinically proven to effectively off-load a minimum of 25% weight-bearing pressure from the forefoot to promote enhanced healing.

The BeneFoot Wedge shoe provides a lower degree of dorsiflexion for enhanced patient comfort whilst continuously reducing forefoot pressure.



BENECARE CODE:

WS

PLEASE TURN TO PAGE 41 FOR SIZING & FITTING OPTIONS

FEATURES & BENEFITS

- Square toe design acts as a bumper to provide additional protection when K-wires are present.
- Advanced tread system increases traction where needed most.
- 15° dorsiflexion to the midfoot and heel.
- Wide base for extra stability.
- Removable forefoot closure provides security whilst eliminating buckle pressure.

INDICATIONS OF USE

- Post-surgical healing.
- Forefoot trauma.
- Wound or ulcerations present on the base of the foot.
- Any situations requiring pressure offload of the forefoot.
- Diabetic foot indications.

BENEKIDZ PAEDIATRIC CAST SHOE

The Benekidz paediatric cast shoe is made from firm, slightly padded upper, with a rigid rocker bottom sole. It benefits from a high-top design to provide added stability around the ankle and compressed foam insole for added comfort.

FEATURES & BENEFITS

- Firm slightly padded upper.
- Rigid Rocker bottom sole.
- 2 wide hook and loop adjustable straps.
- Fits above the ankle for added support and protection.
- Compressed foam insole.

INDICATIONS OF USE

- Over a cast.
- Over a bandage.
- Where normal footwear is prohibited.



BENECARE CODE:

PCS

PLEASE TURN TO PAGE 41 FOR SIZING & FITTING OPTIONS

BENEFOOT B-COMFY SHOE



NEW
PRODUCT
AVAILABLE
SUMMER
2023

BENECARE CODE: BCS

PLEASE TURN TO PAGE 41 FOR SIZING & FITTING OPTIONS

The BeneFoot B-Comfy shoe is suitable for both indoor and outdoor use. The upper is made from a breathable soft material, the whole forefoot of the shoe lifts away and the heel of the shoe opens to allow easy access for people with limited dexterity. The sole of the shoe has a durable nonslip sole and the shoe also benefits from a removable insole that can be replaced for a bespoke insole if clinically required. This lightweight comfy shoe is designed for use where normal footwear is not appropriate e.g diabetic foot ulcers, lymphoedema, and oedematous legs. It will provide support and comfort when mobilising.

FEATURES & BENEFITS

- Soft comfy upper. • Forefoot opening. • Heel opening. • Supportive durable sole. • Removable insole. • Hook and loop closure system.

INDICATIONS OF USE

- Diabetic foot ulcers.
- Lymphoedema.
- Any situation where size of foot prohibits normal footwear.

BENEFOOT B-COMFY BOOT

The BeneFoot B-Comfy Boot is suitable for both indoor and outdoor use. The upper is made from a breathable soft material extended up the ankle to provide added support. The whole forefoot of the boot lifts away and the heel of the boot opens to allow easy access for people with limited dexterity. The sole of the boot has a durable nonslip sole and the boot also benefits from a removable insole that can be replaced for a bespoke insole if clinically required. This lightweight comfy boot is designed for use where normal footwear is not appropriate e.g diabetic foot ulcers, lymphoedema, and oedematous legs. It will provide support and comfort when mobilising.



NEW
PRODUCT
AVAILABLE
SUMMER
2023

FEATURES & BENEFITS

- Soft Comfy upper. • Forefoot opening. • Heel opening. • Supportive durable sole. • Removable insole. • Hook and loop closure system.

INDICATIONS OF USE

- Diabetic foot ulcers.
- lymphoedema.
- Any situation where size of foot prohibits normal footwear.

BENECARE CODE: BCB

PLEASE TURN TO PAGE 41 FOR SIZING & FITTING OPTIONS

BENEFOOT ECO WALKER



BENECARE CODE: FW

PLEASE TURN TO PAGE 41 FOR SIZING & FITTING OPTIONS

The BeneFoot Eco Walker Boot is an ultra-lightweight walker boot providing support and comfort to patients requiring short term immobilisation. The Eco Walker Boot comprises of a rocker sole and removable foam lining, helping to improve short term rehabilitation capabilities

FEATURES & BENEFITS

- Hook & Loop closure provides easy adjustment.
- Ultra low profile design eliminates height difference and facilitates natural gait.
- Fixed stirrups immobilise the ankle securely and comfortably.
- Rocker bottom sole allows easier ambulation.
- Padded soft foam keeps lower leg comfortable.
- Ultra lightweight design makes it suitable for short term immobilisation.

INDICATIONS OF USE

- Replaces the need for a heavy cast.
- Stable fracture of fibula, calcaneus and malleolus.
- Post ligament, soft tissue and tendon surgery.
- Post-surgery management.
- Metatarsals & tarsals.

BENEFOOT FIXED WALKER

The BeneFoot Fixed Walker Boot is strong, durable and benefits from a deluxe terry towelling liner, providing extra comfort for the limb. Strong metal struts with a shock absorbing rocker bottom sole, facilitates a more natural walking gait. It replaces the need for multiple cast changes.

A lightweight Fixed Walker boot provides support and comfort to patients.



FEATURES & BENEFITS

- Hook & Loop closure provides easy adjustment.
- Ultra low profile design eliminates height difference and facilitates natural gait.
- Fixed stirrups immobilise the ankle securely and comfortably.
- Impact resistant and shock absorbing sole for enhanced protection.
- Rocker Bottom sole allows easier ambulation.
- Padded soft foam keeps lower leg comfortable.
- Excellent durability

INDICATIONS OF USE

- Post ligament, soft tissue & tendon surgery.
- Replaces the need for heavy casts.
- Stable fracture of fibula, calcaneus and malleolus.
- Metatarsals & tarsals.
- Post-surgery management.

BENECARE CODE:

FXPRO

PLEASE TURN TO PAGE 41 FOR SIZING & FITTING OPTIONS

BENEFOOT ROM WALKER

The BeneFoot Rom Walker, is a high quality walker boot with a range of motion hinge at the ankle, which benefits from a soft foam liner with strong rigid metal struts. The rocker bottom sole facilitates a more natural walking gait, whilst the shock absorbing sole provides enhanced protection.



BENECARE CODE:

FXROM

PLEASE TURN TO PAGE 42 FOR SIZING & FITTING OPTIONS

FEATURES & BENEFITS

- Replaces the need for multiple cast changes.
- Rom hinge with 10° increments at ankle.
- Soft foam liner for added comfort.
- Rocker bottom sole facilitate more natural walking gait.
- Shock absorbing sole.
- Strong metal struts provide stability of the ankle.
- Hook and loop adjustable fastening system.

INDICATIONS OF USE

- Achilles Tendon Rupture.
- Fracture management of the foot & ankle.
- Post ligament and tendon surgery.

BENEKIDZ WALKER

The Benekidz Paediatric Walker is a lightweight walker boot providing support and comfort for children. The hook and loop closure provides easy adjustment and the fixed stirrups immobilises the ankle securely and comfortably. The rocker bottom sole not only allows for easy ambulation, but is also impact resistant and shock absorbing, with increased protection for the patient. The boot is extremely durable and its ultra-low profile design eliminates height difference and facilitates natural gait. The padded soft foam also ensures the lower leg is constantly comfortable.

FEATURES & BENEFITS

- Rigid metal frame provides protection for active children.
- Shock absorbing sole to improve comfort & protects the user from slips.
- Rocker bottom sole allows easier ambulation.
- Lightweight design for kids to help them keep active.
- Hook & Loop closure provides easy adjustment.
- Ultra low profile design eliminates height difference and facilitates natural gait.
- Padded soft foam keeps lower leg comfortable.
- Excellent durability.

INDICATIONS OF USE

- Stable fracture of fibula, calcaneus & malleolus.
- Metatarsals & tarsals.
- Soft tissue, tendon & ligament injuries.
- Post-surgery management.
- Toe walking.



BENECARE CODE:

BKW

PLEASE TURN TO PAGE 42 FOR SIZING & FITTING OPTIONS

BENEFOOT SWEDISH AFO

The BeneFoot Swedish AFO has a comfortable, low profile prefabricated leaf spring design to correct flaccid foot drop, especially on post-stroke patients. Injection moulded polypropylene allows for different thicknesses throughout the orthosis. Thicker on the vertical aspect for rigidity, thinner footplate for ease of trimming, comfort and flexibility.

**BACK
IN STOCK
SOON**

*Summer 2023

FEATURES & BENEFITS

- Injection-moulded for a high quality finish.
- Easily trimmed footplate.
- May be heat modified for minor fitting adjustments.

INDICATIONS OF USE

- Drop foot.
- Mild drop foot secondary to other neurological pathologies.



BENECARE CODE:

SAFO

PLEASE TURN TO PAGE 42 FOR SIZING & FITTING OPTIONS

BENEFOOT FPFS SUPPORT



The BeneFoot FPFS reduces pain allowing the patient to do strengthening exercises to the foot and ankle in particular the spring ligament of the foot, therefore preventing further discomfort from this condition.

FEATURES & BENEFITS

- Simple and inexpensive way to reduce inflammation.
- Elastic tension helps relieve pain & promote muscle development.
- Anatomically designed to provide comfort and control.
- Available in two adjustable sizes.
- Maintains tension continuously during wear, allowing the patient to exercise.
- Fast, easy 15-second application.
- Breathable stiff foam, wickable fabric draws perspiration away from skin.

INDICATIONS OF USE

- Plantar Fasciitis.
- Inflammation of the spring ligament.
- Painful heel.

BENECARE CODE:

FPFS

PLEASE TURN TO PAGE 42 FOR SIZING & FITTING OPTIONS

ACHILLES HEEL PADS

The BeneFoot Achilles Heel Pad provides additional Achilles tendon protection.

The heel supports can be placed on top of each other, inside a walker boot, to achieve the desired heel raise. Three heel supports equal 22° of heel raise. Two heel supports equal 16° of heel raise. One heel support equals 10° of heel raise.

FEATURES & BENEFITS

- Adjustable height.
- Achieve desired heel raise.
- Minimal patient discomfort.

INDICATIONS OF USE

- Achilles Tendon rupture.
- Tendon and ligament injuries.



BENECARE CODE:

FXSW

PLEASE TURN TO PAGE 42 FOR SIZING & FITTING OPTIONS

B Mobile

TELESCOPIC ROM BRACES

The B Mobile product range is an exciting and innovative range of telescopic knee and elbow rom braces. This dynamic range will assist in the rehabilitation process and is suitable for a variety of indications, including trauma and orthopaedic, post-operative and sporting injuries.



The Lightweight structure assists with the activities of daily living.

10 degree incremental hinge will provide controlled range of motion at the joint.

Extendable struts ensures excellent fit and support to the limb.

Quick release fastening system aids easy application.

Rubberised foam pads provide comfort and prevent migration of the brace.



PRODUCT INFORMATION

BENECARE CODE	PRODUCT NAME	LEFT OR RIGHT	GO TO....
EABL	EXTENDER ARM BRACE	LEFT	Page 25 for more information.
EABR	EXTENDER ARM BRACE	RIGHT	Page 25 for more information.
EABDL	EXTENDER ARM BRACE DELUXE	LEFT	Page 25 for more information.
EABDR	EXTENDER ARM BRACE DELUXE	RIGHT	Page 25 for more information.
EABDGL	EXTENDER ARM BRACE DELUXE (WITH GRIP)	LEFT	Page 25 for more information.
EABDGR	EXTENDER ARM BRACE DELUXE (WITH GRIP)	RIGHT	Page 25 for more information.
EKBD	EXTENDER KNEE BRACE DELUXE	UNIVERSAL	Page 21 for more information.

The B-Mobiille range wiillll help rehabiillitate the patient to a maxiimum attainable degree..

ANKLE BRACES



A sprain of the ankle is a tear to one or more ligaments of the ankle, unfortunately it can be quite difficult to identify the severity and get an accurate diagnosis. This type of injury can often be managed successfully with rest, elevation and a supportive brace.

The BeneCare range of ankle braces provide excellent stability for this type of indication, the ideal choice to promote rehabilitation, they are also an excellent choice to provide support during sporting activities.

EXO ANKLE BRACE

19

MEMORY ANKLE BRACE

19



EXO ANKLE BRACE



The BeneCare Exo Ankle Brace technology forms a plastic exoskeleton directly onto a soft breathable fabric, combining compression with comfort and support. The brace is designed to provide moderate stability and fits easily into a sports shoe so the ideal choice for sports injuries.

FEATURES & BENEFITS

- Plastic exoskeleton provides moderate stability.
- Breathable fabric for comfort.
- Lace strapping allows a custom fit.
- Will fit into sports shoe.

INDICATIONS OF USE

- Inversion/eversion injuries.
- Acute sports injuries.
- Ankle strains and sprains.



BENECARE CODE:

EBAB

PLEASE TURN TO PAGE 42 FOR SIZING & FITTING OPTIONS

MEMORY ANKLE BRACE

The BeneCare Memory Ankle Brace is low profile, designed with a strong plastic shell to provide excellent stability and lightweight memory foam for extra comfort. The soft adjustable strapping system allows a custom fit circumferentially and across the heel providing comfort and support.

**BACK
IN STOCK
SOON**

*Summer 2023

FEATURES & BENEFITS

- Lightweight strong plastic shell for excellent stability.
- Breathable memory foam for added comfort
- Soft adjustable strap with Y tab fasteners.
- Adjustable heel strap.
- Should be worn inside firm shoe or trainer.

INDICATIONS OF USE

- Inversion and Eversion injuries.
- Collateral ligament injuries.
- Mild to moderate ankle sprains.
- Post-operative indications.
- Post cast stabilisation.
- Acute sports injuries.



BENECARE CODE:

MAB

PLEASE TURN TO PAGE 42 FOR SIZING & FITTING OPTIONS

KNEE SUPPORTS



The knee joint is one of the strongest and most important joints in the human body. The complexity of the joint demands a brace to provide stability appropriate to the joint and the indication.

The Benecare range of knee braces are the ideal choice to enable the patient to stay mobile whilst having the comfort and stability required to rehabilitate the injury to the maximum attainable degree.

EXTENDER KNEE BRACE DELUXE	21
-----------------------------------	-----------

HINGED KNEE STABILISER	21
-------------------------------	-----------

ROM KNEE BRACE	22
-----------------------	-----------

TRI PANEL KNEE SPLINT	22
------------------------------	-----------



NEO WRAP KNEE BRACE (SHORT)	23
------------------------------------	-----------

NEO WRAP KNEE BRACE (LONG)	23
-----------------------------------	-----------

PRO ACTION KNEE BRACE	23
------------------------------	-----------



B-MOBILE EXTENDER KNEE BRACE DELUXE

The BeneCare B-Mobile Extender Knee Brace Deluxe is sleek, low profile and functional, as well as helping to promote rehabilitation for a variety of trauma and orthopaedic conditions.

The telescopic struts provide excellent stability, whilst the ROM hinge allows controlled flexion or extension which has 10° increments ranging from 0° to 120°, and the 4 adjustable cuffs can be positioned independently will ensure a custom fit. There are quick release fastening buckles for ease of use and the rubberised non-slip pads prevent migration of the brace and maintain stability.



FEATURES & BENEFITS

- Low profile brace.
- Telescopic struts.
- Quick release fastening system.
- ROM with 10° increments.
- Non-slip pads on the cuffs.
- All cuffs move independently.

INDICATIONS OF USE

- Sports injury.
- Soft tissue injuries of the knee.
- ACL/PCL.
- Post - Operative support.
- Fracture management.

BENECARE CODE: EKBD

PLEASE TURN TO PAGE 42 FOR SIZING & FITTING OPTIONS

HINGED KNEE STABILISER

The BeneCare Hinged Knee Stabiliser benefits from a blend of soft neoprene that will provide soothing warmth by retaining the body heat, whilst the removable polycentric hinge provides stability and prevents hyperextension of the knee. The patella buttress will stabilise and prevent subluxation of the patella.



FEATURES & BENEFITS

- Polycentric hinge will provide medial and lateral stability & also prevent hyperextension.
- Soft blend of neoprene provides soothing warmth by retaining the body heat.
- Opposing adjustable straps to the back and front will provide controlled compression.
- Contoured construction for comfort and support patella buttress will stabilise the patella.

INDICATIONS OF USE

- Rheumatoid arthritis.
- Mild to moderate hyperextension.
- Subluxation of patella.

BENECARE CODE: HKS

PLEASE TURN TO PAGE 42 FOR SIZING & FITTING OPTIONS

THE BENECARE ROM KNEE BRACE

The BeneCare Rom Knee Brace provides comfort and support whilst providing excellent stability. The soft foam pads provide comfort, and the metal struts provide stability of the knee preventing hyperextension, valgus and varus deviation.

The range of motion hinge can control flexion and extension whilst providing excellent stability to the knee.



FEATURES & BENEFITS

- Soft laminated foam body of the brace.
- A variety of lengths available 16", 20" & 24".
- Extension and immobilisation limitations ranging from 0° to 75°.
- Flexion limitations 0° to 120°.

INDICATIONS OF USE

- Soft tissue injuries of the knee.
- ACL/PCL injuries.
- Collateral ligament injuries.
- Post-operative support.
- Fracture management.

BENECARE CODE:

RKB

PLEASE TURN TO PAGE 42 FOR SIZING & FITTING OPTIONS

THE BENECARE TRI PANEL KNEE SPLINT

The BeneCare Tri Panel Knee Splint – is the ideal choice for all trauma and orthopaedic indications where immobilisation of the knee is required. The soft foam three-panel designed provides a functional and custom fit. The mouldable rigid struts in each panel provide the stability, and the soft foam provides extra comfort.



FEATURES & BENEFITS

- Three piece design with soft foam for providing extra comfort and precise fit.
- Mouldable medial lateral and posterior stays provide stability.
- Hook and loop fastening system.

INDICATIONS OF USE

- Fracture management.
- Soft tissue injuries.
- Patella dislocations.
- Acute sporting injuries.
- Pre & Post-operative immobilisation.

BENECARE CODE:

TPKS

PLEASE TURN TO PAGE 42 FOR SIZING & FITTING OPTIONS

NEO-WRAP KNEE BRACE (SHORT)

The BeneCare Neo-Wrap Knee Brace (short) is a soft neoprene front opening brace, with hyper stop hinges providing medial and lateral stability whilst allowing flexion.

The easy open tabs on the hook and loop straps make the brace the ideal solution for arthritic patients with limited dexterity.

Available in Black & Beige.



FEATURES & BENEFITS

- Front opening soft neoprene.
- Universal.
- Open popliteal area for added comfort when flexing of the knee.
- Removable hyper stop hinge.
- Contrasting proximal and distal strap to prevent rotation.
- Easy grip tabs for patients with poor dexterity.

INDICATIONS OF USE

- Soft tissue injuries of the knee.
- Rheumatoid arthritis.
- Mild to moderate hyperextension.
- Mild to moderate valgus/varus deformity.

BENECARE CODE:

SKKB

PLEASE TURN TO PAGE 42 FOR SIZING & FITTING OPTIONS

NEO-WRAP KNEE BRACE (LONG)

The BeneCare Neo-Wrap Knee Brace (long) is a soft neoprene front opening brace, with long hyper stop hinges providing medial and lateral stability whilst allowing flexion.

The easy open tabs on the hook and loop straps make the brace the ideal solution for arthritic patients with limited dexterity.

Available in Black & Beige.



BENECARE CODE:

LSKKB

PLEASE TURN TO PAGE 43 FOR SIZING & FITTING OPTIONS

FEATURES & BENEFITS

- Front opening soft Neoprene.
- Universal.
- Open popliteal area for added comfort when flexing of the knee.
- Removable hyper stop hinge.
- Contrasting proximal and distal strap to prevent rotation.
- Easy grip tabs for patients with poor dexterity.

INDICATIONS OF USE

- Soft tissue injuries of the knee.
- Rheumatoid arthritis.
- Mild to moderate hyperextension.
- Mild to moderate valgus/varus deformity.

PRO ACTION KNEE BRACE

The BeneCare Pro-Action Knee Brace is an advanced medical device used to increase the knee stability due to an ACL/PCL deficiency and to prevent hyperextension via an adjustable ROM hinge. It is also recommended for low to medium impact sports.

FEATURES & BENEFITS

- 4-leverage design for optimal stabilisation.
- ROM hinges with plug-in flex/ext stops.
- Extension limitation at: 0°, 10°, 20°, 30°.
- Flexion limitation at: 0°, 30°, 45°, 60°, 90°.
- Immobilisation limitation at: 0°, 30°.
- Hypoallergenic silicone-grip strap pads for anti-migration.

INDICATIONS OF USE

- Knee instability due to anterior and/or posterior cruciate ligament.
- (ACL/PCL) deficient.
- To restrict the range of motion after reconstructive surgery.



BENECARE CODE:

PAKB

PLEASE TURN TO PAGE 43 FOR SIZING & FITTING OPTIONS

ARM & SHOULDER



The upper limb includes two complex joints: the shoulder and the elbow. Both joints have very different structures and range of motion due to the type of joints.

The BeneCare upper limb supports range from shoulder immobilisers to telescopic extender rom elbow braces. The shoulder immobiliser can prevent any movement at the shoulder joint, and the extender elbow brace can restrict the range of motion whilst maintaining function of a major joint. They are designed to provide comfort and support for trauma and orthopaedic indications.

EXTENDER ARM BRACE DELUXE 25

EXTENDER ARM DELUXE (WITH GRIP) 25

EXTENDER ARM BRACE 25

POLY ARM SLING 26

SLING ON A ROLL 26



B-MOBILE EXTENDER ARM BRACE DELUXE

The BeneCare B-Mobile Extender Arm Brace Deluxe is a versatile telescopic brace. It benefits from strong metal struts with four flexible cuffs that can be moulded to the contours of the arm, each cuff can be positioned independently for a more custom fit, and benefitting from quick release buckle. A range of motion hinge with 10° increments will promote accurate control of the joint.

FEATURES & BENEFITS

- Strong rigid telescopic struts.
- ROM Hinge in 10° increments.
- 4 flexible cuffs.
- Each cuff can be moved independently for a more custom fit.
- Quick release buckles for ease of application.

INDICATIONS OF USE

- Soft Tissue and Fracture management.
- Post-operative.
- Sports injuries.
- Neurological indications.



BENECARE CODE: EABD

PLEASE TURN TO PAGE 43 FOR SIZING & FITTING OPTIONS



BENECARE CODE: EABDG

PLEASE TURN TO PAGE 43 FOR SIZING & FITTING OPTIONS

B-MOBILE EXTENDER ARM BRACE DELUXE (WITH GRIP)

The BeneCare B-Mobile Extender Arm Brace Deluxe with hand grip is a versatile telescopic brace. It benefits from strong metal struts with four flexible cuffs that can be moulded to the contours of the arm. Each cuff can be positioned independently for a more custom fit with quick release buckles and a soft foam hand grip. A range of motion hinge with 10° increments will promote accurate control of the joint.

FEATURES & BENEFITS

- Strong rigid telescopic struts.
- ROM hinge in 10° increments.
- 4 flexible cuffs.
- Each cuff can be moved independently for a more custom fit.
- Quick release buckles for ease of application.
- Soft foam covered metal hand grip.

INDICATIONS OF USE

- Soft tissue and fracture management.
- Post-operative.
- Sports injuries.
- Neurological indications.

B-MOBILE EXTENDER ARM BRACE

The BeneCare Extender Arm Brace benefits from strong telescopic struts with four flexible cuffs to provide a more custom fit. The quick release buckles assist with ease of application and the range of motion hinge has increments of 10° to ensure accurate control at the joint.

FEATURES & BENEFITS

- Strong rigid telescopic struts.
- ROM hinge in 10° increments.
- 4 flexible cuffs.
- Quick release buckles for ease of application.

INDICATIONS OF USE

- Soft tissue and fracture management.
- Post-operative.
- Sports injuries.



BENECARE CODE: EAB

PLEASE TURN TO PAGE 43 FOR SIZING & FITTING OPTIONS

POLY ARM SLING



The BeneCare Poly Arm Sling is designed from soft foam for maximum comfort, the body of the sling has adjustable straps. It is ideal as a simple sling to support the arm or immobilisation of the shoulder by adding the additional body strap.

FEATURES & BENEFITS

- Soft foam throughout the sling.
- Two adjustable straps.
- Body of the sling is adjustable.
- Can be used as a simple sling or shoulder immobiliser.
- Easy to apply with touch close fastening system.
- Washable.

INDICATIONS OF USE

- Fracture management to the arm or shoulder.
- Support for upper limb cast.
- Fractured clavicle.
- Post-operative upper limb and shoulder indications.

BENECARE CODE:

PAS

PLEASE TURN TO PAGE 43 FOR SIZING & FITTING OPTIONS

SLING ON A ROLL

The BeneSling is a stretch resistant soft foam immobilisation system, ideal where support and stabilisation of the upper limb is required. The non-stretch foam can be cut to length and will hold the limb in the appropriate position for the indication by using the Y-Tab fastening system.

FEATURES & BENEFITS

- Stretch resistant soft foam provides excellent comfort and support.
- 30 Y-Tab hook fasteners allow easy positioning.
- Can be cut longitudinally to fit paediatrics.
- Washable.

INDICATIONS OF USE


- Fracture and soft tissue injuries of the upper limb and shoulder.
- Post-operative management.

BENECARE CODE:

BSR

PLEASE TURN TO PAGE 43 FOR SIZING & FITTING OPTIONS





BENEKIDZ

BeneKidz products are here to provide support, care and comfort for all the little ones. Whether they've fallen off their bike or pulled a muscle having a kickabout, there's a product for all little limb injuries.

This range of products are designed to fit comfortably and offer the utmost protection and care.

Whether their recovery journey is long or short, we have something to support all little limbs and keep the future generation fit and healthy.



'Protecting
Little People'

Product Name	Code	Sizes
Benekidz Post-Op Medical Shoe	POPD	Go To P33
Benekidz Post-Op High Top Medical Shoe	POHTPD	Go To P33
Benekidz Walker Boot - Small	BKWS	Go To P33
Benekidz Walker Boot - Medium	BKWM	Go To P33
Benekidz Walker Boot - Large	BKWL	Go To P33
Benekidz Knee Tri-Panel - XX Small	TPKS10	Go To P34
Benekidz Knee Tri-Panel - X Small	TPKS12	Go To P34
Benekidz Knee Tri-Panel - Small	TPKS14	Go To P34
Benekidz 7 Inch Wrist Support - Beige	NWS7BE-XS	Go To P34
Benekidz 7 Inch Wrist Support - Black	BNWS7-XS	Go To P34
Benekidz 9 Inch Wrist Support - Beige	NWS9BE-XS	Go To P34
Benekidz 9 Inch Wrist Support - Black	BNWS9-XS	Go To P34
Benekidz Wrist/Thumb Support Closed - Beige	NT WSC-XS	Go To P35
Benekidz Wrist/Thumb Support Closed - Black	BNT WSC-XS	Go To P35
Benekidz Wrist/Thumb Support Open - Beige	NT WSW-XS	Go To P35
Benekidz Wrist/Thumb Support Open - Black	BNT WSW-XS	Go To P35
Benekidz Wrist/Thumb Wrap - Beige	NW T W-XS	Go To P35
Benekidz Wrist/Thumb Wrap - Black	BNW T W-XS	Go To P35
Benekidz Abducted Wrist/Thumb - Beige	ATW-XS	Go To P36
Benekidz Abducted Wrist/Thumb - Black	BATW-XS	Go To P36
Benekidz CMC Comfort Splint	CMC-XS	Go To P36
Benekidz Paediatric Cast Shoe	PCS	Go To P36

WRIST & THUMB



The hand and wrist have a complex structure, that is required to undertake intricate actions with precise movements. BeneCare have a versatile range of hand and wrist splints, which are designed to meet the challenging lifestyles of today, whilst allowing functionality and providing relevant comfort and support for the patient, to continue with the activities of daily living. Suitable for trauma and orthopaedic indications, including rheumatoid and osteoarthritis.

PRO UNIVERSAL WRIST SPLINT 29

CMC SPLINT 29

NEO WRAP WRIST SPLINT (7 INCH) 29

NEO WRAP WRIST SPLINT (9 INCH) 30

MALLET FINGER SPLINTS 30



NEO WRAP WRIST/THUMB (OPEN) 30

NEO WRAP WRIST/THUMB (CLOSED) 31

NEO WRAP WRIST/THUMB WRAP 31

ABDUCTED WRIST/THUMB WRAP 31



PRO UNIVERSAL WRIST SPLINT

The BeneCare Pro Universal Wrist Splint is constructed using lightweight and durable fabrics with a solid palmer bar. The splint's unique design allows it to be used on both the left and right hand by re-aligning the palmer bar. The elasticated wrist strap provides additional support and stabilisation of the wrist.

FEATURES & BENEFITS

- Universal Left and Right.
- Lightweight durable construction.
- Removable adjustable metal stay.
- Breathable Lycra lining.
- Contoured for the palmer crease of the hand.
- Wide elasticated strap to provide additional support of the wrist.

INDICATIONS OF USE

- Soft tissue and fracture management.
- Sports injuries.
- Arthritis and tendonitis.



BENECARE CODE:

UWS

PLEASE TURN TO PAGE 43 FOR SIZING & FITTING OPTIONS

CMC SPLINT



The BeneCare Comfort CMC Splint constructed of perforated neoprene and benefits from a terry towelling lining that helps to keep the skin cool and comfortable.

It is designed with a contoured thumb strap that will provide gentle compression to the CMC joint of the thumb. It will help relieve discomfort associated to arthritis, tendonitis and prevent subluxation of the CMC joint.

Please do not tumble dry.

FEATURES & BENEFITS

- Perforated neoprene allows breathability.
- Terry towelling liner.
- Soft neoprene with contoured thumb strap will provide comfort & gentle compression.
- Washable 40°.

INDICATIONS OF USE

- Arthritis at the CMC joint.
- Soft tissue indications of the CMC joint.

BENECARE CODE:

CMS

PLEASE TURN TO PAGE 43 FOR SIZING & FITTING OPTIONS

NEOWRAP WRIST SPLINT (7 INCH)

The BeneCare NeoWrap 7" (18cm) Wrist splint is constructed of soft neoprene, that will provide comfort and support in addition to retaining heat to soothe painful joints. The neoprene can be trimmed for a tailored fit and the metal stay can be moulded to the contours of the hand for a more custom fit.

Available in Beige or Black. Please do not tumble dry.

FEATURES & BENEFITS

- Neoprene allows gentle compression whilst the metal stay provides stability to the wrist.
- Easy to modify the neoprene & the metal stay.
- Easy to apply with the hook and loop fastening system.
- Maintains full function of the fingers.
- Washable 40°.

INDICATIONS OF USE

- Soft tissue and Fracture Management.
- Repetitive sprain injury.
- Rheumatoid arthritis.
- Osteoarthritis.



BENECARE CODE:

NWS7BE/BNWS7

PLEASE TURN TO PAGE 43 FOR SIZING & FITTING OPTIONS



BENECARE CODE: NWS9BE/BNWS9

PLEASE TURN TO PAGE 44 FOR SIZING & FITTING OPTIONS

NEOWRAP WRIST SPLINT (9 INCH)

The BeneCare NeoWrap 9" (23cm) Wrist Splint is constructed of soft neoprene, that will provide comfort and support in addition to retaining heat to soothe painful joints. The neoprene can be trimmed for a tailored fit and the metal stay can be moulded to the contours of the hand for a more custom fit.

Available in Beige or Black . Please do not tumble dry.

FEATURES & BENEFITS

- Neoprene allows gentle compression whilst the metal stay provides stability to the wrist.
- Easy to modify the neoprene & the metal stay.
- Easy to apply with the hook and loop fastening system.
- Maintains full function of the fingers.
- Washable 40°.

INDICATIONS OF USE

- Soft tissue and fracture management.
- Repetitive sprain injury.
- Rheumatoid arthritis.
- Osteoarthritis.

MALLET FINGER SPLINTS

BeneCare Mallet Finger Splints are designed to support the distal joint of finger in extension while permitting unrestricted movement of proximal interphalangeal joint. The BeneCare Mallet Finger Splint, is normally used without padding and is taped to middle phalanx, maintaining the distal phalanx in extension.



FEATURES & BENEFITS

- Made from high quality rigid plastic.
- Ventilated design for greater comfort.
- 8 Sizes available.

INDICATIONS OF USE

- Mallet deformities of the distal phalanx of the finger.
- Fracture of the distal phalanx.

BENECARE CODE: FMS

PLEASE TURN TO PAGE 44 FOR SIZING & FITTING OPTIONS



BENECARE CODE: NT WSW/BNT WSW

PLEASE TURN TO PAGE 44 FOR SIZING & FITTING OPTIONS

NEOWRAP WRIST/THUMB SPLINT (OPEN)

The BeneCare NeoWrap Wrist/Thumb Splint 7" (18cm), is made of soft neoprene and benefits from 2 metal stays, one on the palmar aspect, and the second along the lateral border to provide support for the thumb. Available in two designs (also see closed version), the open thumb holds the thumb in a neutral position with some flexibility of the DIP joint.

Available in Beige or Black. Please do not tumble dry.

FEATURES & BENEFITS

- Wrist and thumb Wrap around support.
- Excellent support for the CMC and MCP joint of the thumb.
- Made from soft neoprene for extra comfort.
- Two adjustable metal stays for a more custom fit.
- Washable at 40° do not tumble dry.

INDICATIONS OF USE

- Fracture & trauma.
- Healing ligament injuries.
- Rheumatoid / Osteoarthritis.

NEOWRAP WRIST/THUMB SPLINT (CLOSED)

The BeneCare NeoWrap Wrist/Thumb Splint 7" (18cm), is made of soft neoprene and benefits from 2 metal stays one on the palmer aspect the second along the lateral border to provide support for the thumb. Available in two designs (also see open version). The closed design holds the thumb firmly in abduction.

Available in Beige or Black. Please do not tumble dry.



FEATURES & BENEFITS

- Wrist and Thumb Wrap around support.
- Excellent support for the CMC and MCP joint of the thumb.
- Made from soft neoprene for extra comfort.
- Two adjustable metal stays for a more custom fit.
- Washable at 40° do not tumble dry.

INDICATIONS OF USE

- Fracture & trauma.
- Healing ligament injuries.
- Rheumatoid / Osteoarthritis.

BENECARE CODE:

NTWS/BNTWSC

PLEASE TURN TO PAGE 45 FOR SIZING & FITTING OPTIONS



NEOWRAP WRIST/THUMB WRAP

The BeneCare NeoWrap Wrist/Thumb Wrap is made of soft neoprene and is designed to support the CMC joint of the thumb and the wrist by providing gentle compression.

Available in Beige or Black. Please do not tumble dry.

FEATURES & BENEFITS

- Soft neoprene provides gentle compression.
- Thumb extension for additional support.
- Can be trimmed for a more custom fit.
- Washable at 40° do not tumble dry.

INDICATIONS OF USE

- CMC ligament sprains & strains.
- Synovitis.
- Rheumatoid / Osteoarthritis.

BENECARE CODE:

NW T W / BNW T W

PLEASE TURN TO PAGE 45 FOR SIZING & FITTING OPTIONS

ABDUCTED WRIST/THUMB WRAP

The BeneCare Abducted Wrist/Thumb Wrap is made of soft neoprene and benefits from an adjustable abduction splint to hold and support the thumb. It is designed to support the CMC joint of the thumb and the wrist by providing gentle compression.

Available in Beige or Black. Please do not tumble dry.



FEATURES & BENEFITS

- Soft neoprene provides gentle compression.
- Thumb extension for additional support.
- Can be trimmed for a more custom fit.
- Washable at 40° do not tumble dry.

INDICATIONS OF USE

- CMC ligament sprains & strains.
- Synovitis.
- Rheumatoid / Osteoarthritis.

BENECARE CODE:

ATW/BATW

PLEASE TURN TO PAGE 45 FOR SIZING & FITTING OPTIONS

BENEKIDZ

BeneKidz products are here to provide support, care, and comfort for all the little ones. Whether they've fallen off their bike or pulled a muscle having a kickabout, there's a product for all little limb injuries.

This range of products are designed to fit comfortably and offer the utmost protection and care.

Whether their recovery journey is long or short, we have something to support all little limbs and keep the future generation fit and healthy.



BENEKIDZ POST OP MEDICAL SHOE 33

BENEKIDZ HIGH TOP MEDICAL SHOE 33

FIXED WALKER BOOT 33

TRI PANEL KNEE SPLINT 34

NEOWRAP WRIST SPLINT (7") 34



NEOWRAP WRIST SPLINT (9") 34

NEOWRAP WRIST/THUMB (CLOSED) 35

NEOWRAP WRIST/THUMB (OPEN) 35

WNEROIWSTR/ATPH UWMRBIS WT/RTHAP 2395

ABDUCTED WRIST THUMB WRAP 36

CMC SPLINT 36

BENEKIDZ PAEDIATRIC CAST SHOE 36

BENEKIDZ POST OP MEDICAL SHOE

BeneKidz Original Post Op Medical Shoe, benefits from a rounded toe design, it has a strong rigid sole with a slight rocker bottom, in addition to the deep tread that will prevent slipping.

The black nylon mesh body of the shoe has a padded collar, with reinforced heel counter with a dual density insole for added protection and comfort.

This shoe is the ideal choice to keep active little people mobile, with a lower limb injury by protecting the cast or bandage during robust activities.



FEATURES & BENEFITS

- Strong rigid sole with excellent tread.
- Slight rocker bottom.
- Breathable nylon mesh body.
- Padded collar and reinforced heel counter.
- Dual density insole.
- Adjustable hook and loop closure system capabilities to distribute plantar pressures.

INDICATIONS OF USE

- Post-Operative .
- Over large dressings and bandages.
- Diabetic foot conditions including plantar wounds.
- Over a lower limb cast.
- Any situation where the size of the foot prohibits normal footwear.

BENECARE CODE:

POPD

PLEASE TURN TO PAGE 46 FOR SIZING & FITTING OPTIONS

BENEKIDZ HIGH TOP MEDICAL SHOE



Benefits from black nylon mesh body which is extended up the ankle with a reinforced heel counter and additional strap to prevent slippage and provide extra support.

It is designed with a rounded toe, it has a strong rigid sole with a slight rocker bottom, in addition to the deep tread that will prevent slipping and a dual density insole for added protection and comfort.

This shoe is the ideal choice to keep little people mobile, with a lower limb injury by protecting the cast or bandage during robust activities.

FEATURES & BENEFITS

- Strong rigid sole with excellent tread.
- Slight rocker bottom.
- Breathable nylon mesh body.
- Ankle strap eliminates heel slippage and forward foot movement.
- Extended ankle support and reinforced heel counter.
- Dual Density insole.
- Adjustable hook and loop closure system.

INDICATIONS OF USE

- Post operative.
- Over large dressings and bandages.
- Diabetic foot conditions including plantar wounds.
- Over a lower limb cast.
- Any situation where the size of the foot prohibits normal footwear.

BENECARE CODE:

POHTPD

PLEASE TURN TO PAGE 46 FOR SIZING & FITTING OPTIONS

BENEKIDZ WALKER BOOT

The BeneKidz Walker boot is lightweight, it has fixed stirrups and a rocker bottom sole. It provides excellent stabilisation at the ankle whilst enabling little people to mobilise safely.

The soft foam liner and the shock absorbing sole provides additional comfort for the delicate limb.

The boot is extremely durable to deal with robust activities and its low profile design helps to facilitates a natural gait, whilst the hook and loop fastening system allows ease of application and removal.



FEATURES & BENEFITS

- Rigid metal frame provides extra protection.
- Shock absorbing sole to improve comfort & protects the users from slips.
- Rocker bottom sole allows easier ambulation
- Lightweight design to help active kids.
- Hook & Loop closure provides easy adjustment
- low profile design eliminates height difference and facilitates natural gait.
- Padded soft foam keeps lower leg comfortable
- Excellent durability.

INDICATIONS OF USE

- Replaces the need for heavy casts.
- Stable fracture of fibula, calcaneus and malleolus surgery, ligament and tendon injuries.
- Post-surgery management.
- Toe walking.

BENECARE CODE:

BKW

PLEASE TURN TO PAGE 46 FOR SIZING & FITTING OPTIONS



BENECARE CODE: TPKS10/TPKS12/TPKS14

PLEASE TURN TO PAGE 46 FOR SIZING & FITTING OPTIONS

BENEKIDZ TRI-PANEL KNEE SPLINT

Tri Panel Knee Splint – is the ideal choice for all trauma and orthopaedic indications where immobilisation of the knee is required for little people. The adjustable soft foam tri-panel design, provides a functional and custom fit. The mouldable rigid struts in each panel provide stability and the soft foam provides extra comfort for the delicate skin on little people.

FEATURES & BENEFITS

- Three- piece design with soft foam for providing extra comfort and precise fit.
- Mouldable medial lateral and posterior stays provide stability.
- Hook and loop fastening system.

INDICATIONS OF USE

- Fracture Management.
- Soft Tissue injuries.
- Patella dislocations.
- Acute sporting injuries.
- Pre & post - operative immobilisation.

BENEKIDZ 7 INCH WRIST SPLINT

The BeneKidz NeoWrap 7 inch wrist splint is constructed of soft neoprene. It will provide comfort and support for the delicate skin of little people.

The Neoprene can be trimmed for a tailored fit and the metal stay can be moulded to the contours of the hand to ensure a more custom fit.

Please do not tumble dry.



BENECARE CODE: NWS7BE-XS/BNWS7BE-XS

PLEASE TURN TO PAGE 46 FOR SIZING & FITTING OPTIONS

FEATURES & BENEFITS

- Neoprene allows gentle compression whilst the metal stay provides stability to the wrist.
- Easy to modify the neoprene & the metal stay
- Easy to apply with the hook and loop fastening system.
- Maintains full function of the fingers whilst providing support to the wrist.
- Washable 40°.

INDICATIONS OF USE

- Soft tissue and Fracture Management.
- Repetitive sprain injury.
- Arthritis.



BENECARE CODE: NWS9BE-XS/BNWS9BE-XS

PLEASE TURN TO PAGE 46 FOR SIZING & FITTING OPTIONS

BENEKIDZ 9 INCH WRIST SPLINT

The BeneKidz NeoWrap 9 inch wrist splint is constructed of soft neoprene, the 9" brace will provide extra stability and comfort for the delicate skin of little people.

FEATURES & BENEFITS

- Neoprene allows gentle compression whilst the metal stay provides stability to the wrist
- Easy to modify the neoprene & the metal stay
- Easy to apply with the hook and loop fastening system
- Maintains full function of the fingers
- Washable 40°

INDICATIONS OF USE

- Soft tissue and Fracture Management
- Repetitive sprain injury
- Arthritis

BENEKIDZ WRIST/THUMB SPLINT (CLOSED)

The BeneCare NeoWrap Wrist/Thumb Splint (Closed) is made of soft neoprene. It benefits from 2 metal stays one on the palmer aspect, the second along the lateral border to provide comfort and support for the thumb and wrist. Ideal for the robust activities of little people.

It is available in 2 designs, the open thumb holds the thumb in a neutral position with some flexibility of the DIP joint and closed design holds the thumb firmly in abduction.

FEATURES & BENEFITS

- Wrist and thumb Wrap around support.
- Excellent support for the CMC and MCP joint of the thumb.
- Made from soft neoprene for extra comfort.
- Two adjustable metal stays for a more custom fit.
- Washable at 40° do not tumble dry.

INDICATIONS OF USE

- Fracture management to the wrist hand and thumb.
- CMC ligament sprains & strains.
- Synovitis.
- Arthritis.



BENECARE CODE: BNT WSC-XS/NT WSC-XS

PLEASE TURN TO PAGE 46 FOR SIZING & FITTING OPTIONS

BENEKIDZ WRIST/THUMB SPLINT (OPEN)

The BeneCare NeoWrap Wrist/Thumb Splint (Open) is made of soft neoprene. It benefits from 2 metal stays one on the palmer aspect, the second along the lateral border to provide comfort and support for the thumb and wrist. It is ideal for the robust activities of little people.

It is available in 2 designs. The open thumb holds the thumb in a neutral position with some flexibility of the DIP joint, and closed design holds the thumb firmly in abduction.

Please do not tumble dry.



BENECARE CODE: BNT WSW-XS/NT WSW-XS

PLEASE TURN TO PAGE 46 FOR SIZING & FITTING OPTIONS

FEATURES & BENEFITS

- Wrist and thumb Wrap around support.
- Excellent support for the CMC and MCP joint of the thumb.
- Made from soft neoprene for extra comfort.
- Two adjustable metal stays for a more custom fit.
- Washable at 40°.

INDICATIONS OF USE

- Fracture management to the wrist hand and thumb.
- CMC ligament sprains & strains.
- Synovitis.
- Arthritis.

BENEKIDZ WRIST/THUMB WRAP

The BeneCare NeoWrap Wrist/Thumb Wrap is made of soft neoprene, providing extra comfort for the little people and is designed to support the CMC joint of the thumb and the wrist by providing gentle compression.

Please do not tumble dry.

FEATURES & BENEFITS

- Soft Neoprene provides gentle compression.
- Thumb extension for additional support.
- Can be trimmed for a more custom fit.
- Washable at 40° do not tumble dry.

INDICATIONS OF USE

- CMC ligament sprains & strains.
- Synovitis.
- Rheumatoid / Osteoarthritis.



BENECARE CODE: BNW T W-XS/NW T W-XS

PLEASE TURN TO PAGE 46 FOR SIZING & FITTING OPTIONS



BENECARE CODE: ATW-XS/BATW-XS

PLEASE TURN TO PAGE 46 FOR SIZING & FITTING OPTIONS

BENEKIDZ ABDUCTED WRIST/THUMB WRAP

The Benecare NeoWrap Wrist/Thumb Wrap is made of soft neoprene providing comfort for little people. It also benefits from a metal adjustable abduction splint on the lateral border of the thumb. It is designed to provide support to the CMC joint of the thumb and the wrist by providing gentle compression whilst holding the thumb in abduction.

Please do not tumble dry.

FEATURES & BENEFITS

- Soft Neoprene provides gentle compression.
- Thumb extension for additional support.
- Can be trimmed for a more custom fit.
- Washable at 40°.
- Do not tumble dry.

INDICATIONS OF USE

- CMC ligament sprains & strains.
- Synovitis.
- Rheumatoid / Osteoarthritis.

BENEKIDZ CMC SPLINT

The Benecare Comfort CMC Splint is constructed of perforated neoprene and benefits from a terry towelling lining, which helps to keep the skin of little people cool and comfortable. It is designed with a contoured thumb strap that will provide gentle compression to the CMC joint of the thumb, and therefore relieve discomfort associated to arthritis tendonitis and prevent subluxation of the CMC joint.

Please do not tumble dry.



FEATURES & BENEFITS

- Perforated Neoprene allows breathability.
- Terry towelling liner.
- Soft neoprene with contoured thumb strap will provide comfort and gentle compression and support whilst allowing full function of the hand.
- Washable 40°.

INDICATIONS OF USE

- Arthritis at the CMC joint.
- Soft tissue indications of the CMC joint.

BENECARE CODE: CMC-XS

PLEASE TURN TO PAGE 46 FOR SIZING & FITTING OPTIONS



BENECARE CODE: PCS

PLEASE TURN TO PAGE 46 FOR SIZING & FITTING OPTIONS

BENEKIDZ PAEDIATRIC CAST SHOE

The BeneKidz paediatric cast shoe is made from firm, slightly padded upper, with a rigid rocker bottom sole. It benefits from a high - top design to provide added stability around the ankle and compressed foam insole for added comfort.

FEATURES & BENEFITS

- Firm slightly padded upper.
- Rigid rocker-bottom sole.
- 2 wide hook and loop adjustable straps.
- Fits above the ankle for added support and protection.
- Compressed foam insole.

INDICATIONS OF USE

- Over a cast.
- Over a bandage.
- Where normal footwear is prohibited.

THE BENECARE ACADEMY



Committed to your professional development

Free eLearning Portal

Step By Step Video Tutorials

Downloadable Guidelines

Learn In Your Own Time

Be Trained By Professionals

benecaremedical.com/the-benecare-academy

WATCH ON YOUTUBE
www.youtube.com/benecaremedical



Benecare Direct is our e-commerce website that supplies our healthcare range to the general public, medical professionals and organisations.

SIMPLE & SECURE

Ordering online from Benecare couldn't be easier or safer

EXPERT TECHNICAL SUPPORT

Get expert advice and help from our medical specialists

FIRST CLASS CUSTOMER SERVICE

Our dedicated team is on hand to help with all enquiries

QUALITY ASSURANCE

Products designed, developed and manufactured by us



www.benecaredirect.com

PRODUCT INDEX

CODES	PRODUCT NAME	DESCRIPTION	SIZE	COLOUR	FITTING	PAGE #
CCTB2	Benecast Advanced	Conformable Poly Casting Tape	2"/5cm (Length 3.6m)	Blue	Not Applicable	Page 5
CCTB3	Benecast Advanced	Conformable Poly Casting Tape	3"/7.5cm (Length 3.6m)	Blue	Not Applicable	Page 5
CCTB4	Benecast Advanced	Conformable Poly Casting Tape	4"/10cm (Length 3.6m)	Blue	Not Applicable	Page 5
CC TBK2	Benecast Advanced	Conformable Poly Casting Tape	2"/5cm (Length 3.6m)	Black	Not Applicable	Page 5
CC TBK3	Benecast Advanced	Conformable Poly Casting Tape	3"/7.5cm (Length 3.6m)	Black	Not Applicable	Page 5
CC TBK4	Benecast Advanced	Conformable Poly Casting Tape	4"/10cm (Length 3.6m)	Black	Not Applicable	Page 5
CC TPK2	Benecast Advanced	Conformable Poly Casting Tape	2"/5cm (Length 3.6m)	Pink	Not Applicable	Page 5
CC TPK3	Benecast Advanced	Conformable Poly Casting Tape	3"/7.5cm (Length 3.6m)	Pink	Not Applicable	Page 5
CC TPK4	Benecast Advanced	Conformable Poly Casting Tape	4"/10cm (Length 3.6m)	Pink	Not Applicable	Page 5
CCTP2	Benecast Advanced	Conformable Poly Casting Tape	2"/5cm (Length 3.6m)	Purple	Not Applicable	Page 5
CCTP3	Benecast Advanced	Conformable Poly Casting Tape	3"/7.5cm (Length 3.6m)	Purple	Not Applicable	Page 5
CCTP4	Benecast Advanced	Conformable Poly Casting Tape	4"/10cm (Length 3.6m)	Purple	Not Applicable	Page 5
CCTR2	Benecast Advanced	Conformable Poly Casting Tape	2"/5cm (Length 3.6m)	Red	Not Applicable	Page 5
CCTR3	Benecast Advanced	Conformable Poly Casting Tape	3"/7.5cm (Length 3.6m)	Red	Not Applicable	Page 5
CCTR4	Benecast Advanced	Conformable Poly Casting Tape	4"/10cm (Length 3.6m)	Red	Not Applicable	Page 5
CC T W2	Benecast Advanced	Conformable Poly Casting Tape	2"/5cm (Length 3.6m)	White	Not Applicable	Page 5
CC T W3	Benecast Advanced	Conformable Poly Casting Tape	3"/7.5cm (Length 3.6m)	White	Not Applicable	Page 5
CC T W4	Benecast Advanced	Conformable Poly Casting Tape	4"/10cm (Length 3.6m)	White	Not Applicable	Page 5
CC T W5	Benecast Advanced	Conformable Poly Casting Tape	5"/12.5cm (Length 3.6m)	White	Not Applicable	Page 5
CC TG2	Benecast Advanced	Conformable Poly Casting Tape	2"/5cm (Length 3.6m)	Green	Not Applicable	Page 5
CC TG3	Benecast Advanced	Conformable Poly Casting Tape	3"/7.5cm (Length 3.6m)	Green	Not Applicable	Page 5
CC TG4	Benecast Advanced	Conformable Poly Casting Tape	4"/10cm (Length 3.6m)	Green	Not Applicable	Page 5
FC TB2	Benecast FibreTX	Fibreglass Casting Tape	2"/5cm (Length 3.6m)	Blue	Not Applicable	Page 5
FC TB3	Benecast FibreTX	Fibreglass Casting Tape	3"/7.5cm (Length 3.6m)	Blue	Not Applicable	Page 5
FC TB4	Benecast FibreTX	Fibreglass Casting Tape	4"/10cm (Length 3.6m)	Blue	Not Applicable	Page 5
FC TR2	Benecast FibreTX	Fibreglass Casting Tape	2"/5cm (Length 3.6m)	Red	Not Applicable	Page 5
FC TR3	Benecast FibreTX	Fibreglass Casting Tape	3"/7.5cm (Length 3.6m)	Red	Not Applicable	Page 5
FC TR4	Benecast FibreTX	Fibreglass Casting Tape	4"/10cm (Length 3.6m)	Red	Not Applicable	Page 5
FCTW2	Benecast FibreTX	Fibreglass Casting Tape	2"/5cm (Length 3.6m)	White	Not Applicable	Page 5
FCTW3	Benecast FibreTX	Fibreglass Casting Tape	3"/7.5cm (Length 3.6m)	White	Not Applicable	Page 5
FCTW4	Benecast FibreTX	Fibreglass Casting Tape	4"/10cm (Length 3.6m)	White	Not Applicable	Page 5
FCTW5	Benecast FibreTX	Fibreglass Casting Tape	5"/12.5cm (Length 3.6m)	White	Not Applicable	Page 5
FC TPK2	Benecast FibreTX	Fibreglass Casting Tape	2"/5cm (Length 3.6m)	Pink	Not Applicable	Page 5
FC TPK3	Benecast FibreTX	Fibreglass Casting Tape	3"/7.5cm (Length 3.6m)	Pink	Not Applicable	Page 5
FC TPK4	Benecast FibreTX	Fibreglass Casting Tape	4"/10cm (Length 3.6m)	Pink	Not Applicable	Page 5
FCTG2	Benecast FibreTX	Fibreglass Casting Tape	2"/5cm (Length 3.6m)	Green	Not Applicable	Page 5
FCTG3	Benecast FibreTX	Fibreglass Casting Tape	3"/7.5cm (Length 3.6m)	Green	Not Applicable	Page 5
FCTG4	Benecast FibreTX	Fibreglass Casting Tape	4"/10cm (Length 3.6m)	Green	Not Applicable	Page 5
FC TBK2	Benecast FibreTX	Fibreglass Casting Tape	2"/5cm (Length 3.6m)	Black	Not Applicable	Page 5
FC TBK3	Benecast FibreTX	Fibreglass Casting Tape	3"/7.5cm (Length 3.6m)	Black	Not Applicable	Page 5
FCTP2	Benecast FibreTX	Fibreglass Casting Tape	2"/5cm (Length 3.6m)	Purple	Not Applicable	Page 5
FCTP3	Benecast FibreTX	Fibreglass Casting Tape	3"/7.5cm (Length 3.6m)	Purple	Not Applicable	Page 5
FCTP4	Benecast FibreTX	Fibreglass Casting Tape	4"/10cm (Length 3.6m)	Purple	Not Applicable	Page 5
CTB2	Benecast Select	Polyester Casting Tape	2"/5cm (Length 3.6m)	Blue	Not Applicable	Page 6
CTB3	Benecast Select	Polyester Casting Tape	3"/7.5cm (Length 3.6m)	Blue	Not Applicable	Page 6
CTB4	Benecast Select	Polyester Casting Tape	4"/10cm (Length 3.6m)	Blue	Not Applicable	Page 6
CTR2	Benecast Select	Polyester Casting Tape	2"/5cm (Length 3.6m)	Red	Not Applicable	Page 6
CTR3	Benecast Select	Polyester Casting Tape	3"/7.5cm (Length 3.6m)	Red	Not Applicable	Page 6
CTR4	Benecast Select	Polyester Casting Tape	4"/10cm (Length 3.6m)	Red	Not Applicable	Page 6
CTW2	Benecast Select	Polyester Casting Tape	2"/5cm (Length 3.6m)	White	Not Applicable	Page 6
CTW3	Benecast Select	Polyester Casting Tape	3"/7.5cm (Length 3.6m)	White	Not Applicable	Page 6
CTW4	Benecast Select	Polyester Casting Tape	4"/10cm (Length 3.6m)	White	Not Applicable	Page 6
CTW5	Benecast Select	Polyester Casting Tape	5"/12.5cm (Length 3.6m)	White	Not Applicable	Page 6
CTG2	Benecast Select	Polyester Casting Tape	2"/5cm (Length 3.6m)	Green	Not Applicable	Page 6
CTG3	Benecast Select	Polyester Casting Tape	3"/7.5cm (Length 3.6m)	Green	Not Applicable	Page 6
CTG4	Benecast Select	Polyester Casting Tape	4"/10cm (Length 3.6m)	Green	Not Applicable	Page 6
CTP2	Benecast Select	Polyester Casting Tape	2"/5cm (Length 3.6m)	Purple	Not Applicable	Page 6
CTP3	Benecast Select	Polyester Casting Tape	3"/7.5cm (Length 3.6m)	Purple	Not Applicable	Page 6
CTP4	Benecast Select	Polyester Casting Tape	4"/10cm (Length 3.6m)	Purple	Not Applicable	Page 6
C TPK2	Benecast Select	Polyester Casting Tape	2"/5cm (Length 3.6m)	Pink	Not Applicable	Page 6
C TPK3	Benecast Select	Polyester Casting Tape	3"/7.5cm (Length 3.6m)	Pink	Not Applicable	Page 6

CODES	PRODUCT NAME	DESCRIPTION	SIZE	COLOUR	FITTING	PAGE #
C TPK4	Benecast Select	Polyester Casting Tape	4"/10cm (Length 3.6m)	Pink	Not Applicable	Page 6
C TBK2	Benecast Select	Polyester Casting Tape	2"/5cm (Length 3.6m)	Black	Not Applicable	Page 6
C TBK3	Benecast Select	Polyester Casting Tape	3"/7.5cm (Length 3.6m)	Black	Not Applicable	Page 6
C TBK4	Benecast Select	Polyester Casting Tape	4"/10cm (Length 3.6m)	Black	Not Applicable	Page 6
CTO2	Benecast Select	Polyester Casting Tape	2"/5cm (Length 3.6m)	Orange	Not Applicable	Page 6
CTO3	Benecast Select	Polyester Casting Tape	3"/7.5cm (Length 3.6m)	Orange	Not Applicable	Page 6
CTO4	Benecast Select	Polyester Casting Tape	4"/10cm (Length 3.6m)	Orange	Not Applicable	Page 6
CTY2	Benecast Select	Polyester Casting Tape	2"/5cm (Length 3.6m)	Yellow	Not Applicable	Page 6
CTY3	Benecast Select	Polyester Casting Tape	3"/7.5cm (Length 3.6m)	Yellow	Not Applicable	Page 6
CTY4	Benecast Select	Polyester Casting Tape	4"/10cm (Length 3.6m)	Yellow	Not Applicable	Page 6
					-	
C TSRB2	Benecast FLEX	Flexible Casting Tape	2"/5cm (Length 3.6m)	Blue	Not Applicable	Page 6
C TSRB3	Benecast FLEX	Flexible Casting Tape	3"/7.5cm (Length 3.6m)	Blue	Not Applicable	Page 6
C TSRB4	Benecast FLEX	Flexible Casting Tape	4"/10cm (Length 3.6m)	Blue	Not Applicable	Page 6
C TSRP2	Benecast FLEX	Flexible Casting Tape	2"/5cm (Length 3.6m)	Purple	Not Applicable	Page 6
C TSRP3	Benecast FLEX	Flexible Casting Tape	3"/7.5cm (Length 3.6m)	Purple	Not Applicable	Page 6
C TSRP4	Benecast FLEX	Flexible Casting Tape	4"/10cm (Length 3.6m)	Purple	Not Applicable	Page 6
C TSRPK2	Benecast FLEX	Flexible Casting Tape	2"/5cm (Length 3.6m)	Pink	Not Applicable	Page 6
C TSRPK3	Benecast FLEX	Flexible Casting Tape	3"/7.5cm (Length 3.6m)	Pink	Not Applicable	Page 6
C TSRPK4	Benecast FLEX	Flexible Casting Tape	4"/10cm (Length 3.6m)	Pink	Not Applicable	Page 6
C TSR W2	Benecast FLEX	Flexible Casting Tape	2"/5cm (Length 3.6m)	White	Not Applicable	Page 6
C TSR W3	Benecast FLEX	Flexible Casting Tape	3"/7.5cm (Length 3.6m)	White	Not Applicable	Page 6
C TSR W4	Benecast FLEX	Flexible Casting Tape	4"/10cm (Length 3.6m)	White	Not Applicable	Page 6
C TSRR2	Benecast FLEX	Flexible Casting Tape	2"/5cm (Length 3.6m)	Red	Not Applicable	Page 6
C TSRR3	Benecast FLEX	Flexible Casting Tape	3"/7.5cm (Length 3.6m)	Red	Not Applicable	Page 6
C TSRR4	Benecast FLEX	Flexible Casting Tape	4"/10cm (Length 3.6m)	Red	Not Applicable	Page 6
C TSRBK2	Benecast FLEX	Flexible Casting Tape	2"/5cm (Length 3.6m)	Black	Not Applicable	Page 6
C TSRBK3	Benecast FLEX	Flexible Casting Tape	3"/7.5cm (Length 3.6m)	Black	Not Applicable	Page 6
C TSRG2	Benecast FLEX	Flexible Casting Tape	2"/5cm (Length 3.6m)	Green	Not Applicable	Page 6
C TSRG3	Benecast FLEX	Flexible Casting Tape	3"/7.5cm (Length 3.6m)	Green	Not Applicable	Page 6
C TSRY2	Benecast FLEX	Flexible Casting Tape	2"/5cm (Length 3.6m)	Yellow	Not Applicable	Page 6
C TSRY3	Benecast FLEX	Flexible Casting Tape	3"/7.5cm (Length 3.6m)	Yellow	Not Applicable	Page 6
P.O.P.	Benecast POP	Plaster of Paris Rolls B.P.	5cm x 2.7m	White	Not Applicable	Page 7
2	Benecast POP	Plaster of Paris Rolls B.P.	7.5cm x 2.7m	White	Not Applicable	Page 7
P.O.P.	Benecast POP	Plaster of Paris Rolls B.P.	10cm x 2.7m	White	Not Applicable	Page 7
3	Benecast POP	Plaster of Paris Rolls B.P.	15cm x 2.7m	White	Not Applicable	Page 7
P.O.P.	Benecast POP	Plaster of Paris Rolls B.P.	20cm x 2.7m	White	Not Applicable	Page 7
4						
P.O.P.S	Benecast POP	Plaster of Paris Slabs 5 ply B.P.	10cm x 20m	White	Not Applicable	Page 7
4	Benecast POP	Plaster of Paris Slabs 5 ply B.P.	15cm x 20m	White	Not Applicable	Page 7
P.O.P.S	Benecast POP	Plaster of Paris Slabs 5 ply B.P.	20cm x 20m	White	Not Applicable	Page 7
8					-	
BUP2.5	Benecare Undercast Padding	Synthetic Undercast Padding	5cm x 2.7m	White	Not Applicable	Page 7
BUP3	Benecare Undercast Padding	Synthetic Undercast Padding	7.5cm x 2.7m	White	Not Applicable	Page 7
SUP4	Benecare Undercast Padding	Synthetic Undercast Padding	10cm x 2.7m	White	Not Applicable	Page 7
SUP6	Benecare Undercast Padding	Synthetic Undercast Padding	15cm x 2.7m	White	Not Applicable	Page 7
SR2X15	Benecast Splint-On-A-Roll	Synthetic Splint-On-A-Roll	5cm x 4.5m	White	Not Applicable	Page 7
SR3X15	Benecast Splint-On-A-Roll	Synthetic Splint-On-A-Roll	7.5cm x 4.5m	White	Not Applicable	Page 7
SR4X15	Benecast Splint-On-A-Roll	Synthetic Splint-On-A-Roll	10cm x 4.5m	White	Not Applicable	Page 7
SR5X15	Benecast Splint-On-A-Roll	Synthetic Splint-On-A-Roll	12.5cm x 4.5m	White	Not Applicable	Page 7
SR6X15	Benecast Splint-On-A-Roll	Synthetic Splint-On-A-Roll	15cm x 4.5m	White	Not Applicable	Page 7
PCMS001	Benecast Precut Splints	Colles	Colles	White	Not Applicable	Page 8
PCMS002	Benecast Precut Splints	Short Arm	Short Arm	White	Not Applicable	Page 8
PCMS003	Benecast Precut Splints	Volar	Volar	White	Not Applicable	Page 8
PCMS004	Benecast Precut Splints	Thumb Spica	Thumb Spica	White	Not Applicable	Page 8
PCMS005	Benecast Precut Splints	Long Arm Splint	Long Arm Splint	White	Not Applicable	Page 8
PCMS006	Benecast Precut Splints	Long Leg	Long Leg	White	Not Applicable	Page 8
PCMS007	Benecast Precut Splints	Short Leg	Adult - Short Leg	White	Not Applicable	Page 8
PCMS010	Benecast Precut Splints	Short Leg	Child - Short Leg	White	Not Applicable	Page 8

PRODUCT INDEX

CODES	PRODUCT NAME	DESCRIPTION	SIZE	COLOUR	FITTING	PAGE #
SCCB2	Benecast Crepe Bandage	Cotton Crepe Bandage BP	5cm x 4.5m	Beige	Not Applicable	Page 8
SCCB3	Benecast Crepe Bandage	Cotton Crepe Bandage BP	7.5cm x 4.5m	Beige	Not Applicable	Page 8
SCCB4	Benecast Crepe Bandage	Cotton Crepe Bandage BP	10cm x 4.5m	Beige	Not Applicable	Page 8
SCCB6	Benecast Crepe Bandage	Cotton Crepe Bandage BP	15cm x 4.5m	Beige	Not Applicable	Page 8
SN1	Benecast Stockinette	Cotton Ribbed Stockinette	2.5cm x 20m	Beige	Not Applicable	Page 8
SN2	Benecast Stockinette	Cotton Ribbed Stockinette	5cm x 20m	Beige	Not Applicable	Page 8
SN3	Benecast Stockinette	Cotton Ribbed Stockinette	7.5cm x 20m	Beige	Not Applicable	Page 8
SN4	Benecast Stockinette	Cotton Ribbed Stockinette	10cm x 20m	Beige	Not Applicable	Page 8
SN6	Benecast Stockinette	Cotton Ribbed Stockinette	15cm x 20m	Beige	Not Applicable	Page 8
SN8	Benecast Stockinette	Cotton Ribbed Stockinette	20cm x 20m	Beige	Not Applicable	Page 8
NWSS7.5/25	Benecare Swabs	Non Woven Sterile Swabs	7.5 cm (4 ply)	White	Not Applicable	Page 9
NWS10	Benecare Swabs	Non Woven Swabs	10 x 10 cm (4 ply)	White	Not Applicable	Page 9
BPT7	Benecast Pre-Tape	Pre-	7cm x 20m	White	Applicable	Page 9
BPT10	Benecast Pre-Tape	Tape	10cm x 20m	White	Applicable	Page 9
		Pre-				
POPD	BeneFoot Post-Op Medical Shoe	Original Post-Op Medical Shoe	Paediatric	Black	Paediatric - fits left or right	Page 11
POWS	BeneFoot Post-Op Medical Shoe	Original Post-Op Medical Shoe	UK 2.5 - 4.5 EU 35 - 37.5	Black	Women's Small - fits left or right	Page 11
POWM	BeneFoot Post-Op Medical Shoe	Original Post-Op Medical Shoe	UK 5 - 6.5 EU 38 - 39.5	Black	Women's Medium - fits left or right	Page 11
POWL	BeneFoot Post-Op Medical Shoe	Original Post-Op Medical Shoe	UK 7 - 8.5 EU 40 - 43	Black	Women's Large - fits left or right	Page 11
POMS	BeneFoot Post-Op Medical Shoe	Original Post-Op Medical Shoe	UK 5 - 7 EU 37.5 - 40	Black	Men's Small - fits left or right	Page 11
POMM	BeneFoot Post-Op Medical Shoe	Original Post-Op Medical Shoe	UK 7.5 - 9 EU 41 - 43	Black	Men's Medium - fits left or right	Page 11
POML	BeneFoot Post-Op Medical Shoe	Original Post-Op Medical Shoe	UK 9.5 - 11 EU 44 - 46	Black	Men's Large - fits left or right	Page 11
POMXL	BeneFoot Post-Op Medical Shoe	Original Post-Op Medical Shoe	UK 11.5 - 13 EU 46.5 - 48.5	Black	Men's X-Large - fits left or right	Page 11
POTCWS	BeneFoot Post-Op With Toe Cover	Original Post-Op Medical Shoe W/ Toe Cover	UK 2.5 - 4.5 EU 35 - 37.5	Black	Women's Small - fits left or right	Page 11
POTCWM	BeneFoot Post-Op With Toe Cover	Original Post-Op Medical Shoe W/ Toe Cover	UK 5 - 6.5 EU 38 - 39.5	Black	Women's Medium - fits left or right	Page 11
POTCWL	BeneFoot Post-Op With Toe Cover	Original Post-Op Medical Shoe W/ Toe Cover	UK 7 - 8.5 EU 40 - 43	Black	Women's Large - fits left or right	Page 11
POTCMS	BeneFoot Post-Op With Toe Cover	Original Post-Op Medical Shoe W/ Toe Cover	UK 5 - 7 EU 37.5 - 40	Black	Men's Small - fits left or right	Page 11
POTCMM	BeneFoot Post-Op With Toe Cover	Original Post-Op Medical Shoe W/ Toe Cover	UK 7.5 - 9 EU 41 - 43	Black	Men's Medium - fits left or right	Page 11
POTCML	BeneFoot Post-Op With Toe Cover	Original Post-Op Medical Shoe W/ Toe Cover	UK 9.5 - 11 EU 44 - 46	Black	Men's Large - fits left or right	Page 11
POTCMXL	BeneFoot Post-Op With Toe Cover	Original Post-Op Medical Shoe W/ Toe Cover	UK 11.5 - 13 EU 46.5 - 48.5	Black	Men's X-Large - fits left or right	Page 11
POHTPD	BeneFoot High Top Shoe	High Top Medical Shoe	Paediatric	Black	Paediatric - fits left or right	Page 12
POHT WS	BeneFoot High Top Shoe	High Top Medical Shoe	UK 2.5 - 4.5 EU 35 - 37.5	Black	Women's Small - fits left or right	Page 12
POHT WM	BeneFoot High Top Shoe	High Top Medical Shoe	UK 5 - 6.5 EU 38 - 39.5	Black	Women's Medium - fits left or right	Page 12
POHT WL	BeneFoot High Top Shoe	High Top Medical Shoe	UK 7 - 8.5 EU 40 - 43	Black	Women's Large - fits left or right	Page 12
POHTMS	BeneFoot High Top Shoe	High Top Medical Shoe	UK 5 - 7 EU 37.5 - 40	Black	Men's Small - fits left or right	Page 12
POHTMM	BeneFoot High Top Shoe	High Top Medical Shoe	UK 7.5 - 9 EU 41 - 43	Black	Men's Medium - fits left or right	Page 12
POHTML	BeneFoot High Top Shoe	High Top Medical Shoe	UK 9.5 - 11 EU 44 - 46	Black	Men's Large - fits left or right	Page 12
POHTMXL	BeneFoot High Top Shoe	High Top Medical Shoe	UK 11.5 - 13 EU 46.5 - 48.5	Black	Men's X-Large - fits left or right	Page 12
POHTCWS	BeneFoot High Top Shoe With Toe Cover	High Top Medical Shoe With Toe Cover	UK 2.5 - 4.5 EU 35 - 37.5	Black	Women's Small - fits left or right	Page 12
POHTCWM	BeneFoot High Top Shoe With Toe Cover	High Top Medical Shoe With Toe Cover	UK 5 - 6.5 EU 38 - 39.5	Black	Women's Medium - fits left or right	Page 12
POHTCWL	BeneFoot High Top Shoe With Toe Cover	High Top Medical Shoe With Toe Cover	UK 7 - 8.5 EU 40 - 43	Black	Women's Large - fits left or right	Page 12
POHTCMS	BeneFoot High Top Shoe With Toe Cover	High Top Medical Shoe With Toe Cover	UK 5 - 7 EU 37.5 - 40	Black	Men's Small - fits left or right	Page 12
POHTCMM	BeneFoot High Top Shoe With Toe Cover	High Top Medical Shoe With Toe Cover	UK 7.5 - 9 EU 41 - 43	Black	Men's Medium - fits left or right	Page 12
POHTCML	BeneFoot High Top Shoe With Toe Cover	High Top Medical Shoe With Toe Cover	UK 9.5 - 11 EU 44 - 46	Black	Men's Large - fits left or right	Page 12
POHT	BeneFoot High Top Shoe With Toe Cover	High Top Medical Shoe With Toe Cover	UK 11.5 - 13 EU 46.5 - 48.5	Black	Men's X-Large - fits left or right	Page 12
TCMXL						
BMST WS	BeneFoot Square Toe Shoe	Square Toe Medical Shoe	UK 2.5 - 4.5 EU 35 - 37.5	Black	Women's Small - fits left or right	Page 12
BMST WM	BeneFoot Square Toe Shoe	Square Toe Medical Shoe	UK 5 - 6.5 EU 38 - 39.5	Black	Women's Medium - fits left or right	Page 12
BMST WL	BeneFoot Square Toe Shoe	Square Toe Medical Shoe	UK 7 - 8.5 EU 40 - 43	Black	Women's Large - fits left or right	Page 12
BMSTMS	BeneFoot Square Toe Shoe	Square Toe Medical Shoe	UK 5 - 7 EU 37.5 - 40	Black	Men's Small - fits left or right	Page 12
BMSTMM	BeneFoot Square Toe Shoe	Square Toe Medical Shoe	UK 7.5 - 9 EU 41 - 43	Black	Men's Medium - fits left or right	Page 12
BMSTML	BeneFoot Square Toe Shoe	Square Toe Medical Shoe	UK 9.5 - 11 EU 44 - 46	Black	Men's Large - fits left or right	Page 12
BMSTMXL	BeneFoot Square Toe Shoe	Square Toe Medical Shoe	UK 11.5 - 13 EU 46.5 - 48.5	Black	Men's X-Large - fits left or right	Page 12
BMSTTCWS	BeneFoot Square Toe Shoe With Toe Cover	Square Toe Medical Shoe With Toe Cover	UK 2.5 - 4.5 EU 35 - 37.5	Black	Women's Small - fits left or right	Page 13
BMSTTCWM	BeneFoot Square Toe Shoe With Toe Cover	Square Toe Medical Shoe With Toe Cover	UK 5 - 6.5 EU 38 - 39.5	Black	Women's Medium - fits left or right	Page 13

PRODUCT INDEX

CODES	PRODUCT NAME	DESCRIPTION	SIZE	COLOUR	FITTING	PAGE #
BMSTTCWL	BeneFoot Square Toe Shoe With Toe Cover	Square Toe Medical Shoe With Toe Cover	UK 7 - 8.5 EU 40 - 43	Black	Women's Large - fits left or right	Page 13
BMSTTCMS	BeneFoot Square Toe Shoe With Toe Cover	Square Toe Medical Shoe With Toe Cover	UK 5 - 7 EU 37.5 - 40	Black	Men's Small - fits left or right	Page 13
BMSTTCMM	BeneFoot Square Toe Shoe With Toe Cover	Square Toe Medical Shoe With Toe Cover	UK 7.5 - 9 EU 41 - 43	Black	Men's Medium - fits left or right	Page 13
BMSTTCML	BeneFoot Square Toe Shoe With Toe Cover	Square Toe Medical Shoe With Toe Cover	UK 9.5 - 11 EU 44 - 46	Black	Men's Large - fits left or right	Page 13
BMST	BeneFoot Square Toe Shoe With Toe Cover	Square Toe Medical Shoe With Toe Cover	UK 11.5 - 13 EU 46.5 - 48.5	Black	Men's X-Large - fits left or right	Page 13
TCMXL						
WSS	BeneFoot Wedge Shoe	Off-Loading Wedge Shoe	UK W: 5.5 - 8 M: 5.5 - 7.5 (EU 38 - 41.5)	Black	Small - fits left or right	Page 13
WSM	BeneFoot Wedge Shoe	Off-Loading Wedge Shoe	UK W: 8.5 - 11 M: 8 - 9.5 (EU 42 - 44)	Black	Medium - fits left or right	Page 13
WSL	BeneFoot Wedge Shoe	Off-Loading Wedge Shoe	UK W: 11+ M: 10 - 11.5 (EU 44.5 - 46.5)	Black	Large - fits left or right	Page 13
WSXL	BeneFoot Wedge Shoe	Off-Loading Wedge Shoe	UK Mens 12 - 13.5 (EU 47 - 49)	Black	X-Large - fits left or right	Page 13
PCSS	BeneKidz Paediatric Cast Shoe	Paediatric Cast Shoe	L: 13.6cm W: 5.5cm	Blue	Small - fits left or right	Page 13
PCSM	BeneKidz Paediatric Cast Shoe	Paediatric Cast Shoe	L: 16.7cm W: 6.4cm	Blue	Medium - fits left or right	Page 13
PCSL	BeneKidz Paediatric Cast Shoe	Paediatric Cast Shoe	L: 19.8cm W: 7.3cm	Blue	Large - fits left or right	Page 13
BCS36L	BeneFoot B-Comfy Shoe - Left	Front & Rear Opening Medical Shoe	Size 36	Black	Sold as a pair	Page 14
BCS36R	BeneFoot B-Comfy Shoe - Right	Front & Rear Opening Medical Shoe	Size 36	Black	Sold as a pair	Page 14
BCS37L	BeneFoot B-Comfy Shoe - Left	Front & Rear Opening Medical Shoe	Size 37	Black	Sold as a pair	Page 14
BCS37R	BeneFoot B-Comfy Shoe - Right	Front & Rear Opening Medical Shoe	Size 37	Black	Sold as a pair	Page 14
BCS38L	BeneFoot B-Comfy Shoe - Left	Front & Rear Opening Medical Shoe	Size 38	Black	Sold as a pair	Page 14
BCS38R	BeneFoot B-Comfy Shoe - Right	Front & Rear Opening Medical Shoe	Size 38	Black	Sold as a pair	Page 14
BCS39L	BeneFoot B-Comfy Shoe - Left	Front & Rear Opening Medical Shoe	Size 39	Black	Sold as a pair	Page 14
BCS39R	BeneFoot B-Comfy Shoe - Right	Front & Rear Opening Medical Shoe	Size 39	Black	Sold as a pair	Page 14
BCS41L	BeneFoot B-Comfy Shoe - Left	Front & Rear Opening Medical Shoe	Size 41	Black	Sold as a pair	Page 14
BCS41R	BeneFoot B-Comfy Shoe - Right	Front & Rear Opening Medical Shoe	Size 41	Black	Sold as a pair	Page 14
BCS42L	BeneFoot B-Comfy Shoe - Left	Front & Rear Opening Medical Shoe	Size 42	Black	Sold as a pair	Page 14
BCS42R	BeneFoot B-Comfy Shoe - Right	Front & Rear Opening Medical Shoe	Size 42	Black	Sold as a pair	Page 14
BCS43L	BeneFoot B-Comfy Shoe - Left	Front & Rear Opening Medical Shoe	Size 43	Black	Sold as a pair	Page 14
BCS43R	BeneFoot B-Comfy Shoe - Right	Front & Rear Opening Medical Shoe	Size 43	Black	Sold as a pair	Page 14
BCS45L	BeneFoot B-Comfy Shoe - Left	Front & Rear Opening Medical Shoe	Size 45	Black	Sold as a pair	Page 14
BCS45R	BeneFoot B-Comfy Shoe - Right	Front & Rear Opening Medical Shoe	Size 45	Black	Sold as a pair	Page 14
BCS46L	BeneFoot B-Comfy Shoe - Left	Front & Rear Opening Medical Shoe	Size 46	Black	Sold as a pair	Page 14
BCS46R	BeneFoot B-Comfy Shoe - Right	Front & Rear Opening Medical Shoe	Size 46	Black	Sold as a pair	Page 14
BCS36	BeneFoot B-Comfy Boot - Left	Front & Rear Opening Medical Boot	Size 36	Black	Sold as a pair	Page 14
BCS36	BeneFoot B-Comfy Boot - Right	Front & Rear Opening Medical Boot	Size 36	Black	Sold as a pair	Page 14
BCS37	BeneFoot B-Comfy Boot - Left	Front & Rear Opening Medical Boot	Size 37	Black	Sold as a pair	Page 14
BCS37	BeneFoot B-Comfy Boot - Right	Front & Rear Opening Medical Boot	Size 37	Black	Sold as a pair	Page 14
BCS38	BeneFoot B-Comfy Boot - Left	Front & Rear Opening Medical Boot	Size 38	Black	Sold as a pair	Page 14
BCS38	BeneFoot B-Comfy Boot - Right	Front & Rear Opening Medical Boot	Size 38	Black	Sold as a pair	Page 14
BCS39	BeneFoot B-Comfy Boot - Left	Front & Rear Opening Medical Boot	Size 39	Black	Sold as a pair	Page 14
BCS39	BeneFoot B-Comfy Boot - Right	Front & Rear Opening Medical Boot	Size 39	Black	Sold as a pair	Page 14
BCS41	BeneFoot B-Comfy Boot - Left	Front & Rear Opening Medical Boot	Size 41	Black	Sold as a pair	Page 14
BCS41	BeneFoot B-Comfy Boot - Right	Front & Rear Opening Medical Boot	Size 41	Black	Sold as a pair	Page 14
BCS42	BeneFoot B-Comfy Boot - Left	Front & Rear Opening Medical Boot	Size 42	Black	Sold as a pair	Page 14
BCS42	BeneFoot B-Comfy Boot - Right	Front & Rear Opening Medical Boot	Size 42	Black	Sold as a pair	Page 14
BCS43	BeneFoot B-Comfy Boot - Left	Front & Rear Opening Medical Boot	Size 43	Black	Sold as a pair	Page 14
BCS43	BeneFoot B-Comfy Boot - Right	Front & Rear Opening Medical Boot	Size 43	Black	Sold as a pair	Page 14
BCS45	BeneFoot B-Comfy Boot - Left	Front & Rear Opening Medical Boot	Size 45	Black	Sold as a pair	Page 14
BCS45	BeneFoot B-Comfy Boot - Right	Front & Rear Opening Medical Boot	Size 45	Black	Sold as a pair	Page 14
BCS46	BeneFoot B-Comfy Boot - Left	Front & Rear Opening Medical Boot	Size 46	Black	Sold as a pair	Page 14
BCS46	BeneFoot B-Comfy Boot - Right	Front & Rear Opening Medical Boot	Size 46	Black	Sold as a pair	Page 14
FWS	BeneCare Eco Walker	Fixed Walker Eco - Small	UK 1-3.5 EU 33-36	Black	Fits Left & Right	Page 14
FWM	BeneCare Eco Walker	Fixed Walker Eco - Medium	UK 4-7 EU 37-40	Black	Fits Left & Right	Page 14
FWL	BeneCare Eco Walker	Fixed Walker Eco - Large	UK 7.5-9.5 EU 41-44	Black	Fits Left & Right	Page 14
FXPROS	BeneCare Fixed Walker	Fixed Walker - Small	UK 1-3.5 EU 33-36	Black	Fits Left & Right	Page 15
FXPROM	BeneCare Fixed Walker	Fixed Walker - Medium	UK 4-7 EU 37-40	Black	Fits Left & Right	Page 15
FXPROL	BeneCare Fixed Walker	Fixed Walker - Large	UK 7.5-9.5 EU 41-44	Black	Fits Left & Right	Page 15
FXPROXL	BeneCare Fixed Walker	Fixed Walker - X-Large	UK 10-14 EU 45-48	Black	Fits Left & Right	Page 15

PRODUCT INDEX

CODES	PRODUCT NAME	DESCRIPTION	SIZE	COLOUR	FITTING	PAGE #
FXROMS	BeneCare ROM Walker	ROM Walker - Small	UK 1-3.5 EU 33-36	Black	Fits Left & Right	Page 15
FXROMM	BeneCare ROM Walker	ROM Walker - Medium	UK 4-7 EU 37-40	Black	Fits Left & Right	Page 15
FXROML	BeneCare ROM Walker	ROM Walker - Large	UK 7.5-9.5 EU 41-44	Black	Fits Left & Right	Page 15
FXROMXL	BeneCare ROM Walker	ROM Walker - X-Large	UK 10-14 EU 45-48	Black	Fits Left & Right	Page 15
BKWS	BeneKidz Walker	Paediatric Walker Boot	Small	Black	Fits Left & Right	Page 15
BKWM	BeneKidz Walker	Paediatric Walker Boot	Medium	Black	Fits Left & Right	Page 15
BKWL	BeneKidz Walker	Paediatric Walker Boot	Large	Black	Fits Left & Right	Page 15
SAFOLS	BeneCare Swedish AFO	Swedish AFO Left - Small	UK 1-3.5 EU 33-36	White	Fits small left	Page 16
SAFOLM	BeneCare Swedish AFO	Swedish AFO Left - Medium	UK 4-7 EU 37-40	White	Fits medium left	Page 16
SAFOLL	BeneCare Swedish AFO	Swedish AFO Left - Large	UK 7.5-9.5 EU 41-44	White	Fits large left	Page 16
SAFOLXL	BeneCare Swedish AFO	Swedish AFO Left - X-Large	UK 10-14 EU 45-48	White	Fits x-large left	Page 16
SAFORS	BeneCare Swedish AFO	Swedish AFO Right - Small	UK 1-3.5 EU 33-36	White	Fits small right	Page 16
SAFORM	BeneCare Swedish AFO	Swedish AFO Right - Medium	UK 4-7 EU 37-40	White	Fits medium right	Page 16
SAFORL	BeneCare Swedish AFO	Swedish AFO Right - Large	UK 7.5-9.5 EU 41-44	White	Fits large right	Page 16
SAFORXL	BeneCare Swedish AFO	Swedish AFO Right - X-Large	UK 10-14 EU 45-48	White	Fits x-large right	Page 16
FPFS-SL	BeneFoot FPFS Support	FPFS Support Small - Left	UK 3-5 EU 35-38	Black	Approximate Shoe Size	Page 16
FPFS-LL	BeneFoot FPFS Support	FPFS Support Large - Left	UK 6-14 EU 38.5-49	Black	Approximate Shoe Size	Page 16
FPFS-SR	BeneFoot FPFS Support	FPFS Support Small - Right	UK 3-5 EU 35-38	Black	Approximate Shoe Size	Page 16
FPFS-LR	BeneFoot FPFS Support	FPFS Support Large - Right	UK 6-14 EU 38.5-49	Black	Approximate Shoe Size	Page 16
FXSWL	BeneCare Achilles Heel Pads	Achilles Heel Pads	All Sizes - Left	White	Not Applicable	Page 16
FXSWR	BeneCare Achilles Heel Pads	Achilles Heel Pads	All Sizes - Right	White	Not Applicable	Page 16
EBABS	BeneCare Exo Ankle	Lace-Up Ankle Brace	Small	Black	Circumference: 305 - 330mm, fits left or right	Page 19
EBABM	BeneCare Exo Ankle	Lace-Up Ankle Brace	Medium	Black	Circumference: 330 - 355mm, fits left or right	Page 19
EBABL	BeneCare Exo Ankle	Lace-Up Ankle Brace	Large	Black	Circumference: 355 - 380mm, fits left or right	Page 19
EBABXL	BeneCare Exo Ankle	Lace-Up Ankle Brace	X-Large	Black	Circumference: 380 - 405mm, fits left or right	Page 19
MAB	BeneCare Memory Ankle Brace	Memory Foam Ankle Brace - One Size	Universal	Black	Universal, fits left or right	Page 19
EKBD	B-Mobile Extender Knee Deluxe	Extender Knee Brace Deluxe	Universal	Black	Universal, fits left or right	Page 21
HKSS	BeneCare Hinged Knee Stabiliser	Hinged Knee Stabiliser W/patella	Small 33 - 36cm	Black	Fits left or right. Measure circumference at mid patella	Page 21
HKSM	BeneCare Hinged Knee Stabiliser	Hinged Knee Stabiliser W/patella	Medium 37 - 41cm	Black	Fits left or right. Measure circumference at mid patella	Page 21
HKSL	BeneCare Hinged Knee Stabiliser	Hinged Knee Stabiliser W/patella	Large 42 - 45cm	Black	Fits left or right. Measure circumference at mid patella	Page 21
HKSXL	BeneCare Hinged Knee Stabiliser	Hinged Knee Stabiliser W/patella	X-Large 46 - 50cm	Black	Fits left or right. Measure circumference at mid patella	Page 21
HKSXXL	BeneCare Hinged Knee Stabiliser	Hinged Knee Stabiliser W/patella	XX-Large 51-55cm	Black	Fits left or right. Measure circumference at mid patella	Page 21
RKB16	BeneCare ROM Knee Brace	Soft Foam Padded Knee Brace	Small 16"/41cm	Black	Measure approx 10cm above ankle & 10cm from groin	Page 22
RKB20	BeneCare ROM Knee Brace	Soft Foam Padded Knee Brace	Medium 20"/51cm	Black	Measure approx 10cm above ankle & 10cm from groin	Page 22
RKB24	BeneCare ROM Knee Brace	Soft Foam Padded Knee Brace	Large 24"/61cm	Black	Measure approx 10cm above ankle & 10cm from groin	Page 22
TPKS10	BeneCare Tri-Panel Splint	Universal, fits left or right	10"/25cm	Black	Measure approx 10cm above ankle & 10cm from groin	Page 22
TPKS12	BeneCare Tri-Panel Splint	Universal, fits left or right	12"/30cm	Black	Measure approx 10cm above ankle & 10cm from groin	Page 22
TPKS14	BeneCare Tri-Panel Splint	Universal, fits left or right	14"/35cm	Black	Measure approx 10cm above ankle & 10cm from groin	Page 22
TPKS16	BeneCare Tri-Panel Splint	Universal, fits left or right	16"/41cm	Black	Measure approx 10cm above ankle & 10cm from groin	Page 22
TPKS18	BeneCare Tri-Panel Splint	Universal, fits left or right	18"/46cm	Black	Measure approx 10cm above ankle & 10cm from groin	Page 22
TPKS20	BeneCare Tri-Panel Splint	Universal, fits left or right	20"/51cm	Black	Measure approx 10cm above ankle & 10cm from groin	Page 22
TPKS22	BeneCare Tri-Panel Splint	Universal, fits left or right	22"/56cm	Black	Measure approx 10cm above ankle & 10cm from groin	Page 22
TPKS24	BeneCare Tri-Panel Splint	Universal, fits left or right	24"/61cm	Black	Measure approx 10cm above ankle & 10cm from groin	Page 22
SKKBXS	BeneCare NeoWrap (Short) Brace	Knee Support - Short	X-Small 25-28cm	Beige	Fits left or right. Measure circumference of knee	Page 23
SKKBS	BeneCare NeoWrap (Short) Brace	Knee Support - Short	Small 29-34cm	Beige	Fits left or right. Measure circumference of knee	Page 23
SKKBM	BeneCare NeoWrap (Short) Brace	Knee Support - Short	Medium 35-40cm	Beige	Fits left or right. Measure circumference of knee	Page 23
SKKBL	BeneCare NeoWrap (Short) Brace	Knee Support - Short	Large 41-46cm	Beige	Fits left or right. Measure circumference of knee	Page 23
SKKBXL	BeneCare NeoWrap (Short) Brace	Knee Support - Short	X-Large 47-56cm	Beige	Fits left or right. Measure circumference of knee	Page 23
SKKBXXL	BeneCare NeoWrap (Short) Brace	Knee Support - Short	XX-Large 57-61cm	Beige	Fits left or right. Measure circumference of knee	Page 23
BSKKBXS	BeneCare NeoWrap (Short) Brace	Knee Support - Short	X-Small 25-28cm	Black	Fits left or right. Measure circumference of knee	Page 23
BSKKBBS	BeneCare NeoWrap (Short) Brace	Knee Support - Short	Small 29-34cm	Black	Fits left or right. Measure circumference of knee	Page 23
BSKKBMB	BeneCare NeoWrap (Short) Brace	Knee Support - Short	Medium 35-40cm	Black	Fits left or right. Measure circumference of knee	Page 23
BSKKBBL	BeneCare NeoWrap (Short) Brace	Knee Support - Short	Large 41-46cm	Black	Fits left or right. Measure circumference of knee	Page 23

CODES	PRODUCT NAME	DESCRIPTION	SIZE	COLOUR	FITTING	PAGE #
BSKKBXL	BeneCare NeoWrap (Short) Brace	Knee Support - Short	X-Large 47-56cm	Black	Fits left or right. Measure circumference of knee	Page 23
BSKKBXXL	BeneCare NeoWrap (Short) Brace	Knee Support - Short	XX-Large 57-61cm	Black	Fits left or right. Measure circumference of knee	Page 23
LSKKBXS	BeneCare NeoWrap (Long) Brace	Knee Support - Long	X-Small 25-28cm	Beige	Fits left or right. Measure circumference of knee	Page 23
LSKKBXS	BeneCare NeoWrap (Long) Brace	Knee Support - Long	Small 29-34cm	Beige	Fits left or right. Measure circumference of knee	Page 23
LSKKBXM	BeneCare NeoWrap (Long) Brace	Knee Support - Long	Medium 35-40cm	Beige	Fits left or right. Measure circumference of knee	Page 23
LSKKBXL	BeneCare NeoWrap (Long) Brace	Knee Support - Long	Large 41-46cm	Beige	Fits left or right. Measure circumference of knee	Page 23
LSKKBXL	BeneCare NeoWrap (Long) Brace	Knee Support - Long	X-Large 47-56cm	Beige	Fits left or right. Measure circumference of knee	Page 23
LSKKBXXL	BeneCare NeoWrap (Long) Brace	Knee Support - Long	XX-Large 57-61cm	Beige	Fits left or right. Measure circumference of knee	Page 23
BLSKKBXS	BeneCare NeoWrap (Long) Brace	Knee Support - Long	X-Small 25-28cm	Black	Fits left or right. Measure circumference of knee	Page 23
BLSKKBXS	BeneCare NeoWrap (Long) Brace	Knee Support - Long	Small 29-34cm	Black	Fits left or right. Measure circumference of knee	Page 23
BLSKKBXM	BeneCare NeoWrap (Long) Brace	Knee Support - Long	Medium 35-40cm	Black	Fits left or right. Measure circumference of knee	Page 23
BLSKKBXL	BeneCare NeoWrap (Long) Brace	Knee Support - Long	Large 41-46cm	Black	Fits left or right. Measure circumference of knee	Page 23
BLSKKBXL	BeneCare NeoWrap (Long) Brace	Knee Support - Long	X-Large 47-56cm	Black	Fits left or right. Measure circumference of knee	Page 23
BLSKKBXXL	BeneCare NeoWrap (Long) Brace	Knee Support - Long	XX-Large 57-61cm	Black	Fits left or right. Measure circumference of knee	Page 23
PAKBLS	BeneCare Pro-Action Knee Brace	Pro-action Knee Brace Left	Small	Black	Measure circumference of knee 15.5-18"/39-45cm	Page 23
PAKBLM	BeneCare Pro-Action Knee Brace	Pro-action Knee Brace Left	Medium	Black	Measure circumference of knee 18.5-21"/47-53.5cm	Page 23
PAKBLL	BeneCare Pro-Action Knee Brace	Pro-action Knee Brace Left	Large	Black	Measure circumference of knee 21.5-23.5"/54-60cm	Page 23
PAKBLXL	BeneCare Pro-Action Knee Brace	Pro-action Knee Brace Left	X-Large	Black	Measure circumference of knee 24-26.5"/61-67cm	Page 23
PAKBRS	BeneCare Pro-Action Knee Brace	Pro-action Knee Brace Right	Small	Black	Measure circumference of knee 15.5-18"/39-45cm	Page 23
PAKBRM	BeneCare Pro-Action Knee Brace	Pro-action Knee Brace Right	Medium	Black	Measure circumference of knee 18.5-21"/47-53.5cm	Page 23
PAKBRL	BeneCare Pro-Action Knee Brace	Pro-action Knee Brace Right	Large	Black	Measure circumference of knee 21.5-23.5"/54-60cm	Page 23
PAKBRXL	BeneCare Pro-Action Knee Brace	Pro-action Knee Brace Right	X-Large	Black	Measure circumference of knee 24-26.5"/61-67cm	Page 23
EABDL	BeneCare Extender Arm Deluxe	Extender Arm Brace Deluxe Left	Universal	Black	One Size - Adjustable	Page 25
EABDR	BeneCare Extender Arm Deluxe	Extender Arm Brace Deluxe Right	Universal	Black	One Size - Adjustable	Page 25
EABDGL	Extender Arm Deluxe (Grip)	Extender Arm Brace Deluxe Left with Grip	Universal	Black	One Size - Adjustable	Page 25
EABDGR	Extender Arm Deluxe (Grip)	Extender Arm Brace Deluxe Right with Grip	Universal	Black	One Size - Adjustable	Page 25
EABL	BeneCare Extender Arm	Extender Arm Brace Left	Universal	Black	One Size - Adjustable	Page 25
EABR	BeneCare Extender Arm	Extender Arm Brace Right	Universal	Black	One Size - Adjustable	Page 25
PAS	BeneCare Poly Arm Sling	Poly Arm Sling Deluxe	Universal	Black	Fits left or right	Page 26
BSR	BeneSling	Sling-on-a-Roll (2 rolls)	5.5cmx12m	Blue	Adjustable with Y Tab fastener	Page 26
UWSS	BeneCare Pro Universal Wrist	Universal Wrist Splint	Small	Black	Wrist Circumference 14 - 15.5cm	Page 29
USWM	BeneCare Pro Universal Wrist	Universal Wrist Splint	Medium	Black	Wrist Circumference 16 - 18cm	Page 29
UWSL	BeneCare Pro Universal Wrist	Universal Wrist Splint	Large	Black	Wrist Circumference 18.5 - 20cm	Page 29
UWSXL	BeneCare Pro Universal Wrist	Universal Wrist Splint	X-Large	Black	Wrist Circumference 20.5 - 22cm	Page 29
CMC-XSL	BeneCare CMC Splint	CMC Comfort Thumb Splint Left	X-Small	Black	Wrist Circumference 11.5 - 13.5cm	Page 29
CMC-SL	BeneCare CMC Splint	CMC Comfort Thumb Splint Left	Small	Black	Wrist Circumference 14 - 16cm	Page 29
CMC-ML	BeneCare CMC Splint	CMC Comfort Thumb Splint Left	Medium	Black	Wrist Circumference 16.5 - 19.5cm	Page 29
CMC-LL	BeneCare CMC Splint	CMC Comfort Thumb Splint Left	Large X-	Black	Wrist Circumference 20 - 21.5cm	Page 29
CMC-XLL	BeneCare CMC Splint	CMC Comfort Thumb Splint Left	Large	Black	Wrist Circumference 22 - 24cm	Page 29
CMC-XXLL	BeneCare CMC Splint	CMC Comfort Thumb Splint Left	XX-Large	Black	Wrist Circumference 24.5 - 27cm	Page 29
CMC-XSR	BeneCare CMC Splint	CMC Comfort Thumb Splint Right	X-Small	Black	Wrist Circumference 11.5 - 13.5cm	Page 29
CMC-SR	BeneCare CMC Splint	CMC Comfort Thumb Splint Right	Small	Black	Wrist Circumference 14 - 16cm	Page 29
CMC-MR	BeneCare CMC Splint	CMC Comfort Thumb Splint Right	Medium	Black	Wrist Circumference 16.5 - 19.5cm	Page 29
CMC-LR	BeneCare CMC Splint	CMC Comfort Thumb Splint Right	Large X-	Black	Wrist Circumference 20 - 21.5cm	Page 29
CMC-XLR	BeneCare CMC Splint	CMC Comfort Thumb Splint Right	Large	Black	Wrist Circumference 22 - 24cm	Page 29
CMC-XXLR	BeneCare CMC Splint	CMC Comfort Thumb Splint Right	XX-Large	Black	Wrist Circumference 24.5 - 27cm	Page 29
NWS7BE-XSL	NeoWrap Wrist Splint (7 Inch)	Wrist Support 7" Left	X-Small	Beige	Wrist Circumference 11.5 - 13.5cm	Page 29
NWS7BE-SL	NeoWrap Wrist Splint (7 Inch)	Wrist Support 7" Left	Small	Beige	Wrist Circumference 14 - 16cm	Page 29
NWS7BE-ML	NeoWrap Wrist Splint (7 Inch)	Wrist Support 7" Left	Medium	Beige	Wrist Circumference 16.5 - 19.5cm	Page 29
NWS7BE-LL	NeoWrap Wrist Splint (7 Inch)	Wrist Support 7" Left	Large	Beige	Wrist Circumference 20 - 21.5cm	Page 29
NWS7BE-XLL	NeoWrap Wrist Splint (7 Inch)	Wrist Support 7" Left	X-Large	Beige	Wrist Circumference 22 - 24cm	Page 29
NWS7BE-XSR	NeoWrap Wrist Splint (7 Inch)	Wrist Support 7" Right	X-Small	Beige	Wrist Circumference 11.5 - 13.5cm	Page 29

PRODUCT INDEX

CODES	PRODUCT NAME	DESCRIPTION	SIZE	COLOUR	FIT TING	PAGE #
NWS7BE-SR	NeoWrap Wrist Splint (7 Inch)	Wrist Support 7" Right	Small	Beige	Wrist Circumference 14 - 16cm	Page 29
NWS7BE-MR	NeoWrap Wrist Splint (7 Inch)	Wrist Support 7" Right	Medium	Beige	Wrist Circumference 16.5 - 19.5cm	Page 29
NWS7BE-LR	NeoWrap Wrist Splint (7 Inch)	Wrist Support 7" Right	Large	Beige	Wrist Circumference 20 - 21.5cm	Page 29
NWS7BE-XLR	NeoWrap Wrist Splint (7 Inch)	Wrist Support 7" Right	X-Large	Beige	Wrist Circumference 22 - 24cm	Page 29
BNWS7-XSL	NeoWrap Wrist Splint (7 Inch)	Wrist Support 7" Left	X-Small	Black	Wrist Circumference 11.5 - 13.5cm	Page 29
BNWS7-SL	NeoWrap Wrist Splint (7 Inch)	Wrist Support 7" Left	Small	Black	Wrist Circumference 14 - 16cm	Page 29
BNWS7-ML	NeoWrap Wrist Splint (7 Inch)	Wrist Support 7" Left	Medium	Black	Wrist Circumference 16.5 - 19.5cm	Page 29
BNWS7-LL	NeoWrap Wrist Splint (7 Inch)	Wrist Support 7" Left	Large	Black	Wrist Circumference 20 - 21.5cm	Page 29
BNWS7-XLL	NeoWrap Wrist Splint (7 Inch)	Wrist Support 7" Left	X-Large	Black	Wrist Circumference 22 - 24cm	Page 29
BNWS7-XSR	NeoWrap Wrist Splint (7 Inch)	Wrist Support 7" Right	X-Small	Black	Wrist Circumference 11.5 - 13.5cm	Page 29
BNWS7-SR	NeoWrap Wrist Splint (7 Inch)	Wrist Support 7" Right	Small	Black	Wrist Circumference 14 - 16cm	Page 29
BNWS7-MR	NeoWrap Wrist Splint (7 Inch)	Wrist Support 7" Right	Medium	Black	Wrist Circumference 16.5 - 19.5cm	Page 29
BNWS7-LR	NeoWrap Wrist Splint (7 Inch)	Wrist Support 7" Right	Large	Black	Wrist Circumference 20 - 21.5cm	Page 29
BNWS7-XLR	NeoWrap Wrist Splint (7 Inch)	Wrist Support 7" Right	X-Large	Black	Wrist Circumference 22 - 24cm	Page 29
NWS9BE-XSL	NeoWrap Wrist Splint (9 Inch)	Wrist Support 9" Left	X-Small	Beige	Wrist Circumference 11.5 - 13.5cm	Page 30
NWS9BE-SL	NeoWrap Wrist Splint (9 Inch)	Wrist Support 9" Left	Small	Beige	Wrist Circumference 14 - 16cm	Page 30
NWS9BE-ML	NeoWrap Wrist Splint (9 Inch)	Wrist Support 9" Left	Medium	Beige	Wrist Circumference 16.5 - 19.5cm	Page 30
NWS9BE-LL	NeoWrap Wrist Splint (9 Inch)	Wrist Support 9" Left	Large	Beige	Wrist Circumference 20 - 21.5cm	Page 30
NWS9BE-XLL	NeoWrap Wrist Splint (9 Inch)	Wrist Support 9" Left	X-Large	Beige	Wrist Circumference 22 - 24cm	Page 30
NWS9BE-XSR	NeoWrap Wrist Splint (9 Inch)	Wrist Support 9" Right	X-Small	Beige	Wrist Circumference 11.5 - 13.5cm	Page 30
NWS9BE-SR	NeoWrap Wrist Splint (9 Inch)	Wrist Support 9" Right	Small	Beige	Wrist Circumference 14 - 16cm	Page 30
NWS9BE-MR	NeoWrap Wrist Splint (9 Inch)	Wrist Support 9" Right	Medium	Beige	Wrist Circumference 16.5 - 19.5cm	Page 30
NWS9BE-LR	NeoWrap Wrist Splint (9 Inch)	Wrist Support 9" Right	Large	Beige	Wrist Circumference 20 - 21.5cm	Page 30
NWS9BE-XLR	NeoWrap Wrist Splint (9 Inch)	Wrist Support 9" Right	X-Large	Beige	Wrist Circumference 22 - 24cm	Page 30
BNWS9-XSL	NeoWrap Wrist Splint (9 Inch)	Wrist Support 9" Left	X-Small	Black	Wrist Circumference 11.5 - 13.5cm	Page 30
BNWS9-SL	NeoWrap Wrist Splint (9 Inch)	Wrist Support 9" Left	Small	Black	Wrist Circumference 14 - 16cm	Page 30
BNWS9-ML	NeoWrap Wrist Splint (9 Inch)	Wrist Support 9" Left	Medium	Black	Wrist Circumference 16.5 - 19.5cm	Page 30
BNWS9-LL	NeoWrap Wrist Splint (9 Inch)	Wrist Support 9" Left	Large	Black	Wrist Circumference 20 - 21.5cm	Page 30
BNWS9-XLL	NeoWrap Wrist Splint (9 Inch)	Wrist Support 9" Left	X-Large	Black	Wrist Circumference 22 - 24cm	Page 30
BNWS9-XSR	NeoWrap Wrist Splint (9 Inch)	Wrist Support 9" Right	X-Small	Black	Wrist Circumference 11.5 - 13.5cm	Page 30
BNWS9-SR	NeoWrap Wrist Splint (9 Inch)	Wrist Support 9" Right	Small	Black	Wrist Circumference 14 - 16cm	Page 30
BNWS9-MR	NeoWrap Wrist Splint (9 Inch)	Wrist Support 9" Right	Medium	Black	Wrist Circumference 16.5 - 19.5cm	Page 30
BNWS9-LR	NeoWrap Wrist Splint (9 Inch)	Wrist Support 9" Right	Large	Black	Wrist Circumference 20 - 21.5cm	Page 30
BNWS9-XLR	NeoWrap Wrist Splint (9 Inch)	Wrist Support 9" Right	X-Large	Black	Wrist Circumference 22 - 24cm	Page 30
FMS1	BeneCare Mallet Finger Splint	Mallet Finger Splint	Size 1	Beige	DIP Joint Circumference 4.5cm	Page 30
FMS2	BeneCare Mallet Finger Splint	Mallet Finger Splint	Size 2	Beige	DIP Joint Circumference 5cm	Page 30
FMS3	BeneCare Mallet Finger Splint	Mallet Finger Splint	Size 3	Beige	DIP Joint Circumference 5.5cm	Page 30
FMS4	BeneCare Mallet Finger Splint	Mallet Finger Splint	Size 4	Beige	DIP Joint Circumference 5.7cm	Page 30
FMS5	BeneCare Mallet Finger Splint	Mallet Finger Splint	Size 5	Beige	DIP Joint Circumference 6cm	Page 30
FMS5.5	BeneCare Mallet Finger Splint	Mallet Finger Splint	Size 5.5	Beige	DIP Joint Circumference 7cm	Page 30
FMS6	BeneCare Mallet Finger Splint	Mallet Finger Splint	Size 6	Beige	DIP Joint Circumference 7.5cm	Page 30
FMS7	BeneCare Mallet Finger Splint	Mallet Finger Splint	Size 7	Beige	DIP Joint Circumference 8.5cm	Page 30
NT WSW-XSL	NeoWrap Wrist/Thumb Support (Open)	Wrist/Thumb Support (Open) Left	X-Small	Beige	Wrist Circumference 11.5 - 13.5cm	Page 30
NT WSW-SL	NeoWrap Wrist/Thumb Support (Open)	Wrist/Thumb Support (Open) Left	Small	Beige	Wrist Circumference 14 - 16cm	Page 30
NT WSW-ML	NeoWrap Wrist/Thumb Support (Open)	Wrist/Thumb Support (Open) Left	Medium	Beige	Wrist Circumference 16.5 - 19.5cm	Page 30
NT WSW-LL	NeoWrap Wrist/Thumb Support (Open)	Wrist/Thumb Support (Open) Left	Large	Beige	Wrist Circumference 20 - 21.5cm	Page 30
NT WSW-XLL	NeoWrap Wrist/Thumb Support (Open)	Wrist/Thumb Support (Open) Left	X-Large	Beige	Wrist Circumference 22 - 24cm	Page 30
NT WSW-XSR	NeoWrap Wrist/Thumb Support (Open)	Wrist/Thumb Support (Open) Right	X-Small	Beige	Wrist Circumference 11.5 - 13.5cm	Page 30
NT WSW-SR	NeoWrap Wrist/Thumb Support (Open)	Wrist/Thumb Support (Open) Right	Small	Beige	Wrist Circumference 14 - 16cm	Page 30
NT WSW-MR	NeoWrap Wrist/Thumb Support (Open)	Wrist/Thumb Support (Open) Right	Medium	Beige	Wrist Circumference 16.5 - 19.5cm	Page 30
NT WSW-LR	NeoWrap Wrist/Thumb Support (Open)	Wrist/Thumb Support (Open) Right	Large	Beige	Wrist Circumference 20 - 21.5cm	Page 30
NT WSW-XLR	NeoWrap Wrist/Thumb Support (Open)	Wrist/Thumb Support (Open) Right	X-Large	Beige	Wrist Circumference 22 - 24cm	Page 30
BNTWSW-XSL	NeoWrap Wrist/Thumb Support (Open)	Wrist/Thumb Support (Open) Left	X-Small	Black	Wrist Circumference 11.5 - 13.5cm	Page 30
BNT WSW-SL	NeoWrap Wrist/Thumb Support (Open)	Wrist/Thumb Support (Open) Left	Small	Black	Wrist Circumference 14 - 16cm	Page 30
BNT WSW-ML	NeoWrap Wrist/Thumb Support (Open)	Wrist/Thumb Support (Open) Left	Medium	Black	Wrist Circumference 16.5 - 19.5cm	Page 30
BNT WSW-LL	NeoWrap Wrist/Thumb Support (Open)	Wrist/Thumb Support (Open) Left	Large	Black	Wrist Circumference 20 - 21.5cm	Page 30
BNTWSW-XLL	NeoWrap Wrist/Thumb Support (Open)	Wrist/Thumb Support (Open) Left	X-Large	Black	Wrist Circumference 22 - 24cm	Page 30

CODES	PRODUCT NAME	DESCRIPTION	SIZE	COLOUR	FITTING	PAGE #
BNTWSW-XSR	NeoWrap Wrist/Thumb Support (Open)	Wrist/Thumb Support (Open) Right	X-Small	Black	Wrist Circumference 11.5 - 13.5cm	Page 30
BNT WSW-SR	NeoWrap Wrist/Thumb Support (Open)	Wrist/Thumb Support (Open) Right	Small	Black	Wrist Circumference 14 - 16cm	Page 30
BNT WSW-MR	NeoWrap Wrist/Thumb Support (Open)	Wrist/Thumb Support (Open) Right	Medium	Black	Wrist Circumference 16.5 - 19.5cm	Page 30
BNT WSW-LR	NeoWrap Wrist/Thumb Support (Open)	Wrist/Thumb Support (Open) Right	Large	Black	Wrist Circumference 20 - 21.5cm	Page 30
BNTWSW-XLR	NeoWrap Wrist/Thumb Support (Open)	Wrist/Thumb Support (Open) Right	X-Large	Black	Wrist Circumference 22 - 24cm	Page 30
NT WSC-XSL	NeoWrap Wrist/Thumb Support (Closed)	Wrist/Thumb Support (Closed) Left	X-Small	Beige	Wrist Circumference 11.5 - 13.5cm	Page 31
NT WSC-SL	NeoWrap Wrist/Thumb Support (Closed)	Wrist/Thumb Support (Closed) Left	Small	Beige	Wrist Circumference 14 - 16cm	Page 31
NT WSC-ML	NeoWrap Wrist/Thumb Support (Closed)	Wrist/Thumb Support (Closed) Left	Medium	Beige	Wrist Circumference 16.5 - 19.5cm	Page 31
NT WSC-LL	NeoWrap Wrist/Thumb Support (Closed)	Wrist/Thumb Support (Closed) Left	Large	Beige	Wrist Circumference 20 - 21.5cm	Page 31
NT WSC-XLL	NeoWrap Wrist/Thumb Support (Closed)	Wrist/Thumb Support (Closed) Left	X-Large	Beige	Wrist Circumference 22 - 24cm	Page 31
NT WSC-XSR	NeoWrap Wrist/Thumb Support (Closed)	Wrist/Thumb Support (Closed) Right	X-Small	Beige	Wrist Circumference 11.5 - 13.5cm	Page 31
NT WSC-SR	NeoWrap Wrist/Thumb Support (Closed)	Wrist/Thumb Support (Closed) Right	Small	Beige	Wrist Circumference 14 - 16cm	Page 31
NT WSC-MR	NeoWrap Wrist/Thumb Support (Closed)	Wrist/Thumb Support (Closed) Right	Medium	Beige	Wrist Circumference 16.5 - 19.5cm	Page 31
NT WSC-LR	NeoWrap Wrist/Thumb Support (Closed)	Wrist/Thumb Support (Closed) Right	Large	Beige	Wrist Circumference 20 - 21.5cm	Page 31
NT WSC-XLR	NeoWrap Wrist/Thumb Support (Closed)	Wrist/Thumb Support (Closed) Right	X-Large	Beige	Wrist Circumference 22 - 24cm	Page 31
BNTWSC-XSL	NeoWrap Wrist/Thumb Support (Closed)	Wrist/Thumb Support (Closed) Left	X-Small	Black	Wrist Circumference 11.5 - 13.5cm	Page 31
BNT WSC-SL	NeoWrap Wrist/Thumb Support (Closed)	Wrist/Thumb Support (Closed) Left	Small	Black	Wrist Circumference 14 - 16cm	Page 31
BNT WSC-ML	NeoWrap Wrist/Thumb Support (Closed)	Wrist/Thumb Support (Closed) Left	Medium	Black	Wrist Circumference 16.5 - 19.5cm	Page 31
BNT WSC-LL	NeoWrap Wrist/Thumb Support (Closed)	Wrist/Thumb Support (Closed) Left	Large	Black	Wrist Circumference 20 - 21.5cm	Page 31
BNTWSC-XLL	NeoWrap Wrist/Thumb Support (Closed)	Wrist/Thumb Support (Closed) Left	X-Large	Black	Wrist Circumference 22 - 24cm	Page 31
BNTWSC-XSR	NeoWrap Wrist/Thumb Support (Closed)	Wrist/Thumb Support (Closed) Right	X-Small	Black	Wrist Circumference 11.5 - 13.5cm	Page 31
BNT WSC-SR	NeoWrap Wrist/Thumb Support (Closed)	Wrist/Thumb Support (Closed) Right	Small	Black	Wrist Circumference 14 - 16cm	Page 31
BNT WSC-MR	NeoWrap Wrist/Thumb Support (Closed)	Wrist/Thumb Support (Closed) Right	Medium	Black	Wrist Circumference 16.5 - 19.5cm	Page 31
BNT WSC-LR	NeoWrap Wrist/Thumb Support (Closed)	Wrist/Thumb Support (Closed) Right	Large	Black	Wrist Circumference 20 - 21.5cm	Page 31
BNTWSC-XLR	NeoWrap Wrist/Thumb Support (Closed)	Wrist/Thumb Support (Closed) Right	X-Large	Black	Wrist Circumference 22 - 24cm	Page 31
NW T W-XSL	BeneCare NeoWrap	Wrist/Thumb Wrap Left	X-Small	Beige	Wrist Circumference 11.5 - 13.5cm	Page 31
NW T W-SL	BeneCare NeoWrap	Wrist/Thumb Wrap Left	Small	Beige	Wrist Circumference 14 - 16cm	Page 31
NW T W-ML	BeneCare NeoWrap	Wrist/Thumb Wrap Left	Medium	Beige	Wrist Circumference 16.5 - 19.5cm	Page 31
NW T W-LL	BeneCare NeoWrap	Wrist/Thumb Wrap Left	Large	Beige	Wrist Circumference 20 - 21.5cm	Page 31
NW T W-XLL	BeneCare NeoWrap	Wrist/Thumb Wrap Left	X-Large	Beige	Wrist Circumference 22 - 24cm	Page 31
NW T W-XSR	BeneCare NeoWrap	Wrist/Thumb Wrap Right	X-Small	Beige	Wrist Circumference 11.5 - 13.5cm	Page 31
NW T W-SR	BeneCare NeoWrap	Wrist/Thumb Wrap Right	Small	Beige	Wrist Circumference 14 - 16cm	Page 31
NW T W-MR	BeneCare NeoWrap	Wrist/Thumb Wrap Right	Medium	Beige	Wrist Circumference 16.5 - 19.5cm	Page 31
NW T W-LR	BeneCare NeoWrap	Wrist/Thumb Wrap Right	Large	Beige	Wrist Circumference 20 - 21.5cm	Page 31
NW T W-XLR	BeneCare NeoWrap	Wrist/Thumb Wrap Right	X-Large	Beige	Wrist Circumference 22 - 24cm	Page 31
BNWTW-XSL	BeneCare NeoWrap	Wrist/Thumb Wrap Left	X-Small	Black	Wrist Circumference 11.5 - 13.5cm	Page 31
BNW T W-SL	BeneCare NeoWrap	Wrist/Thumb Wrap Left	Small	Black	Wrist Circumference 14 - 16cm	Page 31
BNW T W-ML	BeneCare NeoWrap	Wrist/Thumb Wrap Left	Medium	Black	Wrist Circumference 16.5 - 19.5cm	Page 31
BNW T W-LL	BeneCare NeoWrap	Wrist/Thumb Wrap Left	Large	Black	Wrist Circumference 20 - 21.5cm	Page 31
BNWTW-XLL	BeneCare NeoWrap	Wrist/Thumb Wrap Left	X-Large	Black	Wrist Circumference 22 - 24cm	Page 31
BNWTW-XSR	BeneCare NeoWrap	Wrist/Thumb Wrap Right	X-Small	Black	Wrist Circumference 11.5 - 13.5cm	Page 31
BNW T W-SR	BeneCare NeoWrap	Wrist/Thumb Wrap Right	Small	Black	Wrist Circumference 14 - 16cm	Page 31
BNW T W-MR	BeneCare NeoWrap	Wrist/Thumb Wrap Right	Medium	Black	Wrist Circumference 16.5 - 19.5cm	Page 31
BNW T W-LR	BeneCare NeoWrap	Wrist/Thumb Wrap Right	Large	Black	Wrist Circumference 20 - 21.5cm	Page 31
BNWTW-XLR	BeneCare NeoWrap	Wrist/Thumb Wrap Right	X-Large	Black	Wrist Circumference 22 - 24cm	Page 31
ATW-XSL	BeneCare NeoWrap	Abducted Wrist Thumb Wrap Left	X-Small	Beige	Wrist Circumference 11.5 - 13.5cm	Page 31
ATW-SL	BeneCare NeoWrap	Abducted Wrist Thumb Wrap Left	Small	Beige	Wrist Circumference 14 - 16cm	Page 31
ATW-ML	BeneCare NeoWrap	Abducted Wrist Thumb Wrap Left	Medium	Beige	Wrist Circumference 16.5 - 19.5cm	Page 31
ATW-LL	BeneCare NeoWrap	Abducted Wrist Thumb Wrap Left	Large	Beige	Wrist Circumference 20 - 21.5cm	Page 31
ATW-XLL	BeneCare NeoWrap	Abducted Wrist Thumb Wrap Left	X-Large	Beige	Wrist Circumference 22 - 24cm	Page 31
ATW-XSR	BeneCare NeoWrap	Abducted Wrist Thumb Wrap Right	X-Small	Beige	Wrist Circumference 11.5 - 13.5cm	Page 31
ATW-SR	BeneCare NeoWrap	Abducted Wrist Thumb Wrap Right	Small	Beige	Wrist Circumference 14 - 16cm	Page 31
ATW-MR	BeneCare NeoWrap	Abducted Wrist Thumb Wrap Right	Medium	Beige	Wrist Circumference 16.5 - 19.5cm	Page 31
ATW-LR	BeneCare NeoWrap	Abducted Wrist Thumb Wrap Right	Large	Beige	Wrist Circumference 20 - 21.5cm	Page 31
ATW-XLR	BeneCare NeoWrap	Abducted Wrist Thumb Wrap Right	X-Large	Beige	Wrist Circumference 22 - 24cm	Page 31
BATW-XSL	BeneCare NeoWrap	Abducted Wrist Thumb Wrap Left	X-Small	Black	Wrist Circumference 11.5 - 13.5cm	Page 31
BATW-SL	BeneCare NeoWrap	Abducted Wrist Thumb Wrap Left	Small	Black	Wrist Circumference 14 - 16cm	Page 31
BATW-ML	BeneCare NeoWrap	Abducted Wrist Thumb Wrap Left	Medium	Black	Wrist Circumference 16.5 - 19.5cm	Page 31
BATW-LL	BeneCare NeoWrap	Abducted Wrist Thumb Wrap Left	Large	Black	Wrist Circumference 20 - 21.5cm	Page 31

PRODUCT INDEX

CODES	PRODUCT NAME	DESCRIPTION	SIZE	COLOUR	FITTING	PAGE #
BATW-XLL	BeneCare NeoWrap	Abducted Wrist Thumb Wrap Left	X-Large	Black	Wrist Circumference 22 - 24cm	Page 31
BATW-XSR	BeneCare NeoWrap	Abducted Wrist Thumb Wrap Right	X-Small	Black	Wrist Circumference 11.5 - 13.5cm	Page 31
BATW-SR	BeneCare NeoWrap	Abducted Wrist Thumb Wrap Right	Small	Black	Wrist Circumference 14 - 16cm	Page 31
BATW-MR	BeneCare NeoWrap	Abducted Wrist Thumb Wrap Right	Medium	Black	Wrist Circumference 16.5 - 19.5cm	Page 31
BATW-LR	BeneCare NeoWrap	Abducted Wrist Thumb Wrap Right	Large	Black	Wrist Circumference 20 - 21.5cm	Page 31
BATW-XLR	BeneCare NeoWrap	Abducted Wrist Thumb Wrap Right	X-Large	Black	Wrist Circumference 22 - 24cm	Page 31
POPD	BeneKidz Post-Op Medical Shoe	Paediatric - Post-Op Medical Shoe	12-1/EU 30 - 33.5cm	Black	Fits left or right	Page 33
POHTPD	BeneKidz Post-Op High Top Shoe	Paediatric - Post-Op Medical Shoe	12-1/EU 30 - 33.5cm	Black	Fits left or right	Page 33
BKWS	BeneKidz Walker Boot	Paediatric Walker Boot - Small	UK 5 - 7 EU 22 - 24	Black	Fits left or right	Page 33
BKWM	BeneKidz Walker Boot	Paediatric Walker Boot - Medium	UK 7 - 12 EU 24 - 30	Black	Fits left or right	Page 33
BKWL	BeneKidz Walker Boot	Paediatric Walker Boot - Large	UK 12 - 1 EU 30 - 32.5	Black	Fits left or right	Page 33
TPKS10	BeneKidz Tri-Panel Knee	Paediatric Tri-Panel Knee - XX Small	10" (25.4cm)	Black	Fits left or right	Page 34
TPKS12	BeneKidz Tri-Panel Knee	Paediatric Tri-Panel Knee - X Small	12" (30.5cm)	Black	Fits left or right	Page 34
TPKS14	BeneKidz Tri-Panel Knee	Paediatric Tri-Panel Knee - Small	14" (35.5cm)	Black	Fits left or right	Page 34
NWS7BE-XSL	BeneKidz 7 Inch Wrist Support	Paediatric - 7 Inch Wrist Support - Left	Paediatric	Beige	Wrist Circumference 11.5 - 13.5cm	Page 34
NWS7BE-XSR	BeneKidz 7 Inch Wrist Support	Paediatric - 7 Inch Wrist Support - Right	Paediatric	Beige	Wrist Circumference 11.5 - 13.5cm	Page 34
BNWS7BE-XSL	BeneKidz 7 Inch Wrist Support	Paediatric - 7 Inch Wrist Support - Left	Paediatric	Black	Wrist Circumference 11.5 - 13.5cm	Page 34
BNWS7BE-XSR	BeneKidz 7 Inch Wrist Support	Paediatric - 7 Inch Wrist Support - Right	Paediatric	Black	Wrist Circumference 11.5 - 13.5cm	Page 34
NWS9BE-XSL	BeneKidz 9 Inch Wrist Support	Paediatric - 9 Inch Wrist Support - Left	Paediatric	Beige	Wrist Circumference 11.5 - 13.5cm	Page 34
NWS9BE-XSR	BeneKidz 9 Inch Wrist Support	Paediatric - 9 Inch Wrist Support - Right	Paediatric	Beige	Wrist Circumference 11.5 - 13.5cm	Page 34
BNWS9BE-XSL	BeneKidz 9 Inch Wrist Support	Paediatric - 9 Inch Wrist Support - Left	Paediatric	Black	Wrist Circumference 11.5 - 13.5cm	Page 34
BNWS9BE-XSR	BeneKidz 9 Inch Wrist Support	Paediatric - 9 Inch Wrist Support - Right	Paediatric	Black	Wrist Circumference 11.5 - 13.5cm	Page 34
NT WSC-	BeneKidz Wrist/Thumb Support	Wrist/Thumb Support (Closed) - Left	X-Small	Beige	Wrist Circumference 11.5 - 13.5cm	Page 35
XSL NT	BeneKidz Wrist/Thumb Support	Wrist/Thumb Support (Closed) - Right	X-Small	Beige	Wrist Circumference 11.5 - 13.5cm	Page 35
WSC-XSR						
BNTWSC-XSL	BeneKidz Wrist/Thumb Support	Wrist/Thumb Support (Closed) - Left	X-Small	Black	Wrist Circumference 11.5 - 13.5cm	Page 35
BNTWSC-XSR	BeneKidz Wrist/Thumb Support	Wrist/Thumb Support (Closed) - Right	X-Small	Black	Wrist Circumference 11.5 - 13.5cm	Page 35
NT WSW-	BeneKidz Wrist/Thumb Support	Wrist/Thumb Support (Open) - Left	X-Small	Beige	Wrist Circumference 11.5 - 13.5cm	Page 35
XSL NT	BeneKidz Wrist/Thumb Support	Wrist/Thumb Support (Open) - Right	X-Small	Beige	Wrist Circumference 11.5 - 13.5cm	Page 35
WSW-XSR						
BNTWSW-XSL	BeneKidz Wrist/Thumb Support	Wrist/Thumb Support (Open) - Left	X-Small	Black	Wrist Circumference 11.5 - 13.5cm	Page 35
BNTWSW-XSR	BeneKidz Wrist/Thumb Support	Wrist/Thumb Support (Open) - Right	X-Small	Black	Wrist Circumference 11.5 - 13.5cm	Page 35
NW T W-	BeneKidz Wrist/Thumb Wrap	Paediatric - Wrist/Thumb Wrap - Left	X-Small	Beige	Wrist Circumference 11.5 - 13.5cm	Page 35
XSL NW T	BeneKidz Wrist/Thumb Wrap	Paediatric - Wrist/Thumb Wrap - Right	X-Small	Beige	Wrist Circumference 11.5 - 13.5cm	Page 35
W-XSR						
BNWTW-XSL	BeneKidz Wrist/Thumb Wrap	Paediatric - Wrist/Thumb Wrap - Left	X-Small	Black	Wrist Circumference 11.5 - 13.5cm	Page 35
BNWTW-XSR	BeneKidz Wrist/Thumb Wrap	Paediatric - Wrist/Thumb Wrap - Right	X-Small	Black	Wrist Circumference 11.5 - 13.5cm	Page 35
ATW-XSL	BeneKidz Abucted Wrist/Thumb	Abducted Wrist Thumb Wrap - Left	X-Small	Beige	Wrist Circumference 11.5 - 13.5cm	Page 36
ATW-XSR	BeneKidz Abucted Wrist/Thumb	Abducted Wrist Thumb Wrap - Right	X-Small	Beige	Wrist Circumference 11.5 - 13.5cm	Page 36
BATW-XSL	BeneKidz Abucted Wrist/Thumb	Abducted Wrist Thumb Wrap - Left	X-Small	Black	Wrist Circumference 11.5 - 13.5cm	Page 36
BATW-XSR	BeneKidz Abucted Wrist/Thumb	Abducted Wrist Thumb Wrap - Right	X-Small	Black	Wrist Circumference 11.5 - 13.5cm	Page 36
CMC-XSL	BeneKidz CMC Splint	CMC Comfort Thumb Splint - Left	X-Small	Black	Wrist Circumference 11.5 - 13.5cm	Page 36
CMC-XSR	BeneKidz CMC Splint	CMC Comfort Thumb Splint - Right	X-Small	Black	Wrist Circumference 11.5 - 13.5cm	Page 36
PCSS	BeneKidz Paediatric Cast Shoe	Paediatric Cast Shoe - Small	L: 13.6cm x W: 5.5cm	Blue	Fits left or right	Page 36
PCSM	BeneKidz Paediatric Cast Shoe	Paediatric Cast Shoe - Medium	L: 16.7cm x W: 6.4cm	Blue	Fits left or right	Page 36
PCSL	BeneKidz Paediatric Cast Shoe	Paediatric Cast Shoe - Large	L: 19.8cm x W: 7.3cm	Blue	Fits left or right	Page 36

ECOMED® ANTIMICROBIAL & ANTIVIRAL DISPOSABLE HOSPITAL CURTAINS

SURROUND YOU AND PROTECT YOU.

EcoMed
Disposable Hospital Curtains

PARTNERED WITH
MICROBAN 

Contents

Why Change to EcoMed® Antimicrobial Disposable Curtains?

1

Why Choose Microban® Built-in Antimicrobial Technology?

2

Global Laboratory Testing Reports

4

EcoMed Curtain Instructions

5

EcoMed Curtain Selection Guide

6

Who We Are

9

Why Change to EcoMed® Antimicrobial Disposable Curtains?

Reducing: Healthcare Associated Infections (HAIs)

Traditional fabric curtains, though washed regularly, can retain infectious pathogens, posing a risk of cross-contamination between patients and staff. Disposable curtains, on the other hand, are designed for single use, significantly reducing the chances of infection transmission!

Saving Money


A traditional fabric curtain costs times than a disposable curtain and requires high operational cost in terms of man-hours and quarterly washing, while disposable curtains cost zero on this.

Saving Time

Disposable curtains were much easier to hang or remove due to their lighter weight than traditional fabric curtains.

100% Recyclable

Made of 100% polypropylene non-woven fabric, are 100% recyclable as non-clinical waste to ensure that the disposals are safe and environmentally friendly.



EcoMed antimicrobial disposable curtains are designed to fit most existing hospital tracks, are the quick, economical and hygienic alternative to conventional cubicle curtains.

Reference: Jeremy Lim, Operations Admin and Material Management; Priscilla Han, Infection Control, Nursing; Mark Tan, Operations Admin, National Cancer Centre Singapore. "Risk assessment on use of traditional curtains vs disposable curtains in terms of cost effectiveness, work processes and infection control concerns." Singapore Healthcare Management (2018)

Treated With Microban® Built-in Antimicrobial Technology



**Safe to Use
Long Lasting Protection.**



Unlike disinfectants which provide a limited residual activity once the treated surface dries, built-in Microban antimicrobial technology works to continuously minimize the presence of microbes throughout the entire life-cycle of a curtain.

KEY FACT: Even after disinfection, bacterial numbers rebound by 2 hours!

Integrated Microban protection will remain present in EcoMed disposable curtains for its expected lifetime. The built-in technology will not wash off, wear away or leach out.



Protection for the lifetime of the curtain.

Microban antimicrobial technology is not a spray-on but built-in to the polymer structure of the non-woven fabric at the point of manufacture. The antimicrobial active components are evenly distributed throughout the whole curtain. Its efficacy covers the entire lifespan of the curtain.



Safety of Microban.

Microban technology has undergone extensive independent laboratory testing and has a long history of safe use. The biocidal active components of Microban are notified with the EU Biocidal Products Regulation (BPR), registered with the US Environmental Protection Agency (EPA), US Food and Drug Administration (FDA) and OEKO-TEX listed.



*Microban® is a registered trademark of Microban Products Company.

Microban Ultra-Fresh CA-16 is not formulated with, and does not contain, any of the following substances:

- a) Lead; b) Mercury; c) Cadmium; d) Hexavalent chromium; e) Polybrominated biphenyls (PBBs);
- f) Polybrominated diphenyl ethers (PBDEs); g) Bis(2-ethylhexyl) phthalates (DEHPs); h) Butyl benzyl phthalates (BBPs);
- i) Dibutyl phthalates (DBPs); or j) Diisobutyl phthalates (DIBPs)

How Does Microban® Work? Coating VS Built-in Technology

How Does Microban Work?

Microban® antimicrobial technology works at a cellular level to continually disrupt the growth and reproduction of microorganisms. It operates a multi modal attack, damaging the protein, cell membrane, DNA and internal systems of a microbe. Once built into a curtain, Microban® technology will start to work as soon as a microorganism comes into contact with the protected curtain surface.

How does Microban antimicrobial additive differ from disinfectants?

Unlike disinfectants which provide a limited residual activity once the treated surface dries, Microban® built-in antimicrobial technology works to continuously minimize the presence of microbes throughout the entire life-cycle of a curtain. See the differences between disinfectants coating and Microban® built-in antimicrobial technology applied on disposable curtains in the below table.

	Disinfectants Coating	Microban® Built-in Antimicrobial Technology
Product Form	Liquid Form	Masterbatch
Treating Method	Coating on Surface of Fabric	Built into Polymer Structure of Fabric
Effectiveness	Short & Limited	Long Lasting
Treated Curtain Surface	Wet & Sticky	Dry & Smooth
Wipe Clean	NO	YES
Allergic Risk	YES	NO (Oeko-Tex Listed)
Safe for Children	NO	YES (FDA Approved)

*Disclaimer: The comparative information of disinfectants in the above table are only applicable to disinfectants coated on disposable curtains made of non-woven fabric.

Global Laboratory Testing

Antimicrobial & Antiviral / Fire retardant

Curtain Fabric Testing

EcoMed disposable hospital curtains treated with Microban technology are tested for their antibacterial and antiviral efficacy in accordance with internationally recognised test standards ISO 20743 and ISO18184.

ISO 20743: Determination of antibacterial activity of textile products. ISO 20743 specifies quantitative test methods to determine the antibacterial activity of all antibacterial textile products including non-woven fabric.

ISO 18184: Determination of antiviral activity of textile products. ISO 18184 is an antiviral textile test that measures virucidal antimicrobial activity on textiles, non-woven and other porous materials against specific viruses.

Three Years Testing

Tests for antimicrobial activity have been conducted on EcoMed curtain which was used for three years. The curtain has been proven to still kill up to 99.99% of bacteria.

*Testing reports are available upon request.



Bacteria	International Standards	Log Reduction
E. coli	ISO 20743	5.40
MRSA	ISO 20743	5.60
VRE	ISO 20743	4.30
P. aeruginosa	ISO 20743	5.80
E. faecalis	ISO 20743	4.50
S. aureus	ISO 20743	3.60
A. baumannii	ISO 20743	5.50
L. monocytogenes	ISO 20743	3.20
S. enteritidis	ISO 20743	4.00

Viruses	Antiviral Active Value
Human Coronavirus (HCoV-229E)	99.99%
H1N1 Influenza A virus (Swine Flu)	99.99%
Human enterovirus 71 (EV-A71)	99.99%
Coxsackievirus A21	99.63%
Adenovirus type 3	99.56%

Fire retardant standards

EcoMed disposable curtains have been tested to the following International standards by SGS:

Country	International Standards
United Kingdom	BS 5867 Part 2 Type C: 2008
United States	NFPA 701: 2023
Europe	EN13772, EN13773 Class 1
Australia	AS 2755.2-1985
New Zealand	AS/NZS 1530.3: 1999



Disclaimer: EcoMed disposable curtain has been treated with a Microban® antiviral technology to reduce the viral load on the surface by 99.99% within 2 hours, as per the ISO 18184 test method (against H1N1 Influenza A virus). Testing conducted at SGS. Microban® technology is not designed to protect the user of this product or others against disease-causing microorganisms and is not a substitute for normal cleaning practices. And Microban® antiviral technology is not designed to protect users or others against viruses. Normal cleaning and hygiene practices should be maintained.

EcoMed Curtain Instructions

Recommended changeover Procedure

To get the best from EcoMed disposable curtains, please follow the below steps:

- Clean track - follow manufacturer's instructions
- Lubricate track using silicone wipe
- When not in use, just fold away against the wall using the curtain fastening to avoid soiling.
- Due to built-in antimicrobial technology, you can just wipe the curtain surface to clean dust or stain.
- Curtain should be changed if heavily soiled with blood or body fluid to avoid infection.
- Unless contaminated, then dispose with clinical waste.

Recommended changeover procedure

The following table is a comparative look at EcoMed disposable curtains change over frequency and the recommended changeover of reusable curtains.


Area	Recommended changeover frequency of EcoMed Disposable Curtains	For comparative purposes, recommended changeover of Reusable Curtains
Very High Risk <ul style="list-style-type: none"> • Outbreaks of infectious diseases 	Change immediately and upon discharge of the patient from bed or room and dispose	Change immediately and clean upon discharge
High Risk <ul style="list-style-type: none"> • Intensive care units • High dependency unit • Burns unit, renal units • Day surgery • Operating suites • Emergency wards • Wards for patients with highly transmissible diseases e.g. MROs (VRE, MRSA, S. aureus, E. coli) 	Annually - Change annually or if they are visibly soiled or torn	Monthly - Change or clean monthly or upon discharge if MROs present. If visibly soiled change immediately
Significant Risk <ul style="list-style-type: none"> • General wards 	Annually - Change annually or if they are visibly soiled or torn	Six Monthly - Change or clean twice per year and immediately if visibly soiled
Low Risk <ul style="list-style-type: none"> • Rehabilitation • Long term care/nursing homes • Office based practice • Medical centre, radiography 	Every two years - Change every two years or if they are visibly soiled or torn	Annually - Change or clean annually and immediately if visibly soiled

*Note: EcoMed disposable curtains can be changed more frequently depending upon an individual healthcare facilities requirements or changeover and cleaning policy.

Disposable Curtain Selection Guide

1. Hook Type

Our integrated hooks are designed to fit a range of curtain tracks in any healthcare environment.



2. Colour Range

We offer 9 stocked colours, along with custom pediatric printing patterns and a wide selection of bespoke colours.



3. Curtain Size

Available in four standard widths:
• 7.5m • 6.0m • 4.5m • 3.0m

Available in two standard drops:
• 2.0m for suspended curtain tracks.
• 2.3m including 0.3m mesh on top for ceiling mounted curtain tracks to meet fire codes by allowing light, air and water circulation.

*Bespoke curtain widths and drops can be tailored to your needs on request.



EcoMed Disposable Curtain Specification Sheet



- + EcoMed patented hooks can quick and easy fit onto any existing cubicle tracks.
- + Two permanent heat seal self-auditing labels on each curtain allows anyone to see when the curtain was last changed. And the clear NO ENTRY mark to provide patients privacy!
- + EcoMed disposable curtains are permanently pleated and when not in use can be folded back into a much smaller stack thus showing a neat and modern appearance.
- + Made of 100% recyclable polypropylene non-woven fabric treated with built-in Microban antimicrobial technology, extremely durable and can last for over two years.
- + Flame retardant.
- + Each curtain is fitted with a stitched tieback with VELCRO sticker which secures the curtain neatly when not in use.
- + Once the curtain is hanging, remove the protective, hygienic polythene wrapper.

- + The hook for different curtain track could be different, please contact your local agent for details before ordering.
- + Please visit our website to learn more www.ecomedcurtains.com

EU Declaration of Conformity
NON-STERILIZED DISPOSABLE MEDICAL CURTAINS

Regulation (EU) 2017/745 of the European Parliament and of the Council of 5 April 2017 concerning Medical Devices. The undersigned declares, under their sole responsibility, that the products described in this document meet the Council provisions that apply to them and the CE Mark may be affixed.

General Product Name	NON-STERILIZED DISPOSABLE MEDICAL CURTAINS
Legal Manufacturer:	Changzhou EcoMed Technologies Co., Ltd. Business Operation Address: No. 2 Jianye Road, New District, Changzhou City, Jiangsu Province, P.R. China
Description:	Single use medical curtains packed as bulk non-sterile units.
Intended Purpose:	Disposable medical devices for infection control and prevention.
MDR Classification:	Class 1 as per Rule 1 of Annex VIII
EU Authorised Representative:	Caretechion GmbH Niederrheinstraße. 71, 40474 Düsseldorf, Germany
EU Authorised Representative SRN:	DE-AR-000005946
Medical Device Regulation Assessment Route:	Self-certification by Medical Device Regulation Article 52: Conformity assessment based on technical documentation (Annex II & Annex III), and EU Declaration of Conformity (Annex IV).

Applicable Standards

This present declaration is also in conformity with the following European standards and Common Specifications (CS): 2017/745, ISO 13485:2016+A11:2021, ISO 14971:2019+A11:2021, ISO 15223-1:2021, ISO 20417:2021, BS EN ISO 10993 -1:2020.

Declaration of Conformity

Changzhou EcoMed Technologies Co., Ltd., with sole responsibility for the manufacture, packaging and labelling of the NON-STERILIZED DISPOSABLE MEDICAL CURTAINS listed on this certificate, declares that they conform to the relevant provisions of Regulation (EU) 2017/745 and ISO 13485 registered Quality Management System, regardless of whether these operations are carried out by the manufacturer or on their behalf by a third party.

Date: 18 JANUARY, 2024 **Place:** CHANGZHOU, CHINA **Position:** CEO **Signature**



Who We Are: A Leading Designer & Manufacturer



“At EcoMed, we don’t just offer infection control solutions, we pave the way for a safer, cleaner future in healthcare worldwide.”

Simon Wu, Managing Director

About Us

We started in 2012 with a simple idea that every patient, staff and visitor should have well protection from infection in hospitals or other healthcare facilities. We’re dedicated to applying the latest innovative technology to our disposable curtains to reduce and eliminate the infection unseen because we believe that the hospital curtains with the latest technologies will radically change the way to balance hospital outlay and infection control!

Partnered with Microban, we successfully launched the 3rd generation built-in antimicrobial disposable hospital curtains. Microban built-in antimicrobial technology works to continuously minimize the presence of microbes throughout the entire life-cycle of a curtain. We are proud to see it has unbelievably changed the way to prevent infection control.

Committed to excellence, we strive to lead globally in hospital solutions, prioritizing quality, and customer satisfaction. We meticulously oversee every aspect of manufacturing, from materials to packaging while prioritizing team well-being and environmental sustainability.



PARTNERED WITH
MICROBAN



Wild Goose products are all made of
Genuine Australian sheepskin,
not wool on fabric material.

Medical
Product Range

Medical Product Range

Made from genuine Australian sheepskins, individually selected to have the highest density and finest wool qualities. **Proven to reduce pressure and abrasion on the skin** making it ideal for pressure ulcer treatment and prevention.

These products offer pressure care for those who are immobile, frail or in hospital and are ideal for diabetes patients.

- **Comfort.** Cool in summer and warm in winter and absorbs up to 36% of its dry weight in moisture without feeling damp.
- **Durability.** The HiTemp Medical Sheepskin range is urine & blood resistant
- **Suitable** for use in the home, hospital and nursing home



Elbow Protector MR-301

The Elbow Protector aids in healing and also acts as a useful preventative. Treats and protects a specific problem area.

Colour:

Green

Sizes:

Small, Large



Heel Protector MR-201

The Heel Protector minimises pressure on the heel while leaving most of the foot free. Suitable for those being treated for, or at risk of pressure ulcers.

Colour:

Green

Sizes:

Small, Medium, Large





Wheelchair Cover

MR- 401

Specially designed for wheel chairs, it can also be used as a back rest or seat cover. An ideal choice for regular users, keeping them dry and comfortable. Aids in reducing pressure points.

Colour:

Green

Dimensions:

40W x 93L cm

Padded lining, velcro straps



Day Chair Cover

MR-521

A more luxurious cover with extra length to cover the back and seat of a chair. Made from our medical grade sheepskins and easily transferred between a chair and the car this provide comfort and pressure care for users.

Colour:

Green

Dimensions:

50W x 140L cm

Padded lining, velcro straps



Chair Cover

MR-501

A portable seat comfort can be carried around wherever you travel. An ideal choice for regular users, keeping them dry and comfortable. Aids in reducing pressure points.

Colour:

Green

Dimensions:

40W x 48L cm

Padded lining, velcro straps



Overlay

A comfortable pure Australian medical sheepskin with a material backing to ensure the comfort of those spending long hours in bed. Aids in the prevention and reduction of pressure ulcers by providing even pressure distribution. Keeps the skin dry by offering a high level of moisture adsorption. The user does not suffer irritation caused by perspiration and wetness. Includes velcro strapping to make sure it remains comfortably in place.

Colour:

Green

Dimensions:

MR-125 Single bed 75W x 150L cm
MR-121 Day bed 75W x 125 L cm

Padded lining, velcro straps

Medical Rug / Underlay

MR -101

Suitable for anyone, anywhere. Can be used on the bed, chair or floor. Excellent choice for those spending long hours in bed. Aids in the prevention and reduction of pressure ulcers by providing even pressure distribution. Keeps the skin dry by offering a high level of moisture adsorption. The user does not suffer irritation caused by perspiration and wetness.

Colour:

Green

Dimensions:

60-65 W x 90-95 L cm

Natural shape

Foot Rest Overlay M R -131

A luxuriously padded medical sheepskin foot rest overlay to place on the end of your bed, providing a cost effective solution for pressure care.

Colour:

Green

Dimensions:

30W x 70L cm

Padded lining

Available
in **Mens &
Womens**
Sizing.



Wrap Around Boot M R -711

Adjustable Velcro opening makes the boot easy to put on and take off. The special design allows feet rest in snug and cosy boot.

Colour:

Green

Sizes:

XS, S, M, L, XL

Care Instructions

Washable to 80°C and is warranted to 50 washes in hot water. Wash with woolwash only. Ensure velcro is completely closed to stop wool fluff from sticking to the strap.

Do not dry in direct sunlight, tumble dry on warm setting only

Soft Soled Cosy Boot MR-721

Popular with diabetes sufferers and those requiring pressure care. This soft soled boot with adjustable velcro opening provided with genuine Australian medical sheepskin uppers. Perfect to slip on at home or to provide additional comfort while travelling.

Colour:

Green

Sizes:

XS, S, M, L, XL

Open Toed Care Boot M R -761

An open toed boot with extra adjustability for those recovering from surgery or with swelling. Made from medical grade sheepskins with velcro opening at the top of the foot for easy adjustability to provide for your comfort.

Colour:

Green

Sizes:

XS, S, M, L, XL

Size chart	X Small	Small	Medium	L arge	X Large
AUST Men	3-4	5-6	7-8 9-	9-10	11-12
AUST Women	5-6	7-8	10 7-8	11-12	-
UK	3-4	5-6	8-9 9-	9-10	11-12
US Men	4-5	6-7	10 41-	10 -11	12 -13
US Women	5-6	7-8	42	11-12	-
Europe	37-38	39-40		43-44	45-46

Practitioner Guide
including NHS
Supply Chain
Order Codes

Limbo®

The original waterproof protector

Helping patients shower and bathe with confidence

Waterproof Protection for Casts,
Air Walker Boots, Dressings
and PICC Lines/Midlines



Help restore patients' well-being with a LimbO waterproof protector

Thank you for your professional interest in the LimbO Waterproof Protector. This guide explains why the LimbO is recognised as a leading waterproof protection solution for patients and quickly assists you with product selection. The guide highlights how to obtain a LimbO including options on FP10 prescription, NHS Supply Chain (NHS SC) and the NHS Shared Business Services framework (NHS SBS). We have also included a section about our outdoor cast protectors, cast accessories, mobility aids and a range of diabetic socks.

For further support and resources please visit www.limboproducts.co.uk or see our contact details on the back cover.



Why is the LimbO the most recommended protector?

Now on
NHS Supply
Chain

Helping the Patient

- More independence during recovery
- Easy to use
- Maintain personal hygiene
- Durable and submersible
- Reliable waterproof protection
- Adult and child sizes
- Easy to order via phone, web or Amazon
- Prices from £12.00
- Rapid delivery
- Latex free
- Intended for single patient use

Helping Your Team

- Reduce healthcare risks associated with wet dressings and wet casts
- Reduce avoidable NHS costs resulting from re-casting and re-dressing
- Easy size selection
- Available on FP10 prescription, NHS Supply Chain and NHS SBS framework
- Bespoke sizes available for unique cases
- 29 years as a trusted solution
- Dedicated practitioner content on our website

Contents

Benefits of the LimbO	3
Adult LimbOs	4 - 7
Elbow LimbOs (adult & child)	5
Child LimbOs	8 - 9
Cast / Dressing / PICC Line Accessories	10
Outdoor Cast / Dressing Protectors	11
Daily Living Aids and Diabetic Socks	12
Body Block	13
Size Guide	14 - 15
Contact Details	Back cover

"Patient use of the LimbO reduces the chance of wet dressings and casts and the associated problems of discomfort, clinical complications and delays to recovery".

Featuring the unique LimbO neoprene seal

- Flexible and comfortable
- Superior design providing a complete waterproof barrier against the skin
- Gentle non-constrictive seal, made with premium neoprene
- Guaranteed for twelve weeks with regular use
- Proudly designed and manufactured in the UK



"Since using the LimbO sleeves for our PICC patients we have had fewer PICCs falling out and our infection numbers have fallen."
Meinir Hughes, Velindre NHS Trust -
Intravenous Access Nurse Specialist

Guaranteed waterproof protection for casts, dressings or PICC lines

Adult Half Arm

For casts or dressings below the elbow.

- Seals just above the elbow
- Suitable for the bath or shower
- NHS Supply Chain options



Height (cm)	Weight (kg)	Length (cm)	Circumference (cm) Limb	Cast	Model	NHS SBS	NHS SC Order Code
152+	Under 51	49	22-25	33	M50A		N/A
152+	51-70	53	25-29	40	M60	✓	EAZ 706
152+	70+	53	30-39	44	M67	✓	EAZ 699

Adult Full Arm

For casts or dressings covering the whole arm.

- Seals at the top of the upper arm
- Suitable for the bath or shower
- NHS Supply Chain options



Height (cm)	Weight (kg)	Length (cm)	Circumference (cm) Limb	Cast	Model	NHS SBS	NHS SC Order Code
152+	51-70	69	25-29	40	M60L	✓	EAZ 716
152+	70-108	69	30-39	44	M70	✓	EAZ 715
152+	108+	69	39-54	57	M87	✓	EAZ 703

Adult and Child Elbow

For PICC Lines, Midlines, or dressings on the elbow or upper arm.

- Seals above and below the elbow
- Ideal to support patients during IV treatment
- Suitable for the bath or shower
- Prescription options
- Child option
- NHS Supply Chain options



Height (cm)	Weight (kg)	Length (cm)	Circumference (cm) Limb	Dressing	Model	PIP Code	NHS SBS	NHS SC Order Code
Under 152	Under 51	12	22-25	33	M45S	N/A		N/A
152+	Under 51	26	22-25	33	M45	N/A		N/A
152+	51-70	26	25-29	40	M65	386-4303	✓	EAZ 705
152+	70-108	26	30-39	44	M75	386-4311	✓	EAZ 711
152+	108+	26	39-54	57	M85	N/A		N/A

Adult Hand

For light dressings on the hand.

- Seals on the forearm
- Not suitable for plaster casts or following thumb surgery. Select Adult Half Arm model instead



Height (cm)	Weight (kg)	Length (cm)	Circumference (cm) Limb	Cast	Model
Slim to Medium Build		30	22-25	Bandage Only	MITT

The comfortable, trusted solution for staying dry in the shower or bath



Adult Half Leg

For casts or dressings below the knee.

- Seals just above the knee
- Suitable for use with diabetic ulcers
- Suitable for the bath or shower
- Prescription options
- NHS Supply Chain options



Height (cm)	Weight (kg)	Length (cm)	Circumference (cm) Limb	Cast	Model	PIP Code	NHS SBS	NHS SC Order Code
165 to 183	64-102	73	39-54	57	M80	320-0045	✓	EAZ 718
Under 165	64-102	65	39-54	57	M80S	320-0052	✓	EAZ 707
183+	64-102	83	39-54	57	M80L	N/A		N/A
165 to 183	Under 64	73	35-39	44	M76	333-7581	✓	EAZ 710
Under 165	Under 64	65	35-39	44	M76S	333-7599	✓	EAZ 702
183+	Under 64	83	35-39	44	M76L	N/A		N/A
165 to 183	102-140	73	52-57	59	M180	320-0060	✓	EAZ 709
Under 165	102-140	65	52-57	59	M180S	320-0078	✓	EAZ 717
183+	102-140	83	52-57	59	M180L	N/A		N/A
152 to 183	140+	73	65-90	94	M200HL	N/A	✓	EAZ 719



Adult Full Leg

For casts or dressings covering the whole leg.

- Seals at the top of the leg
- Suitable for use with diabetic ulcers
- Suitable for the bath or shower
- NHS Supply Chain options



Height (cm)	Weight (kg)	Length (cm)	Circumference (cm) Limb	Cast	Model	NHS SBS	NHS SC Order Code
165 to 183	64-102	96	52-65	70	M100	✓	EAZ 720
Under 165	64-102	85	52-65	70	M100S		N/A
183+	64-102	103	52-65	70	M100L		N/A
Under 183	Under 64	96	39-54	57	M80LL	✓	EAZ 722
152 to 183	102+	96	65-90	94	M200FL	✓	EAZ 721



Adult Knee

For light dressings on the knee.

- Seals above and below the knee
- For shower use only



Height (cm)	Weight (kg)	Length (cm)	Circumference (cm) Limb	Dressing	Model
165 to 183	64-102	38	39-54	57	M86

Adult Orthopaedic Boot

For orthopaedic boots. For the seal to be large enough to go over the boot, please select the appropriate size Adult Full Leg LimbO model.

- Seals at the top of the thigh
- Suitable for the bath or shower



Adult Foot

For light dressings below the ankle.

- Seals above the ankle bone
- Suitable for use with diabetic ulcers
- For shower use only
- Available on FP10 prescription
- NHS Supply Chain options



Shoe Size	Circumference (cm) Limb	Dressing	Model	PIP Code	NHS SBS	NHS SC Order Code
Max shoe size 11	20-25	33	M20	416-6476	✓	EAZ 723
Max shoe size 11	23-34	40	M25	416-6484	✓	EAZ 701

"Good quality product, good price and delivery service."
Emma Southwell, Orthopaedic Practitioner,
Queen Elizabeth Hospital, King's Lynn

Eliminate the need for replacing wet bandages or ruined plaster casts



Child Full Leg

For casts or dressings covering the leg.

- Seals at the top of the thigh
- Suitable for the bath or shower
- NHS Supply Chain options



Age	Weight (kg)	Height	Length (cm)	Circumference (cm) Limb	Cast	Model	NHS SC Order Code
3 yrs			53	25-34	40	FL3	N/A
4-5 yrs	Size is typical for age, select as you would for clothing		58	25-34	40	FL45	N/A
6-7 yrs			61	35-39	44	FL67	EAZ 714
8-10 yrs			73	39-54	57	FL810	EAZ 704
11-13 yrs			83	39-54	57	FL1113	EAZ 700



Child Half Leg

For casts or dressings below the knee, on older children.

- Seals just above the knee
- Suitable for the bath or shower



Age	Weight (kg)	Height	Length (cm)	Circumference (cm) Limb	Cast	Model
8-10 yrs	Size is typical for age, select as you would for clothing		58	25-34	40	BK810
11-13 yrs			61	35-39	44	BK1113



Child Full Arm

For casts or dressings on the arm.

- Seals at the top of the arm
- Aqua blue junior style
- Suitable for the bath or shower
- NHS Supply Chain options



Age	Weight (kg)	Height	Length (cm)	Circumference (cm) Limb	Cast	Model	NHS SC Order Code
2-3 yrs			31	16-22	29	FA23	N/A
4-5 yrs	Size is typical for age, select as you would for clothing		43	16-22	29	FA45	N/A
6-7 yrs			49	22-25	33	FA67	EAZ 712
8-10 yrs			58	22-25	33	FA810	EAZ 708
11-13 yrs			58	25-34	40	FA1113	EAZ 713



Child Half Arm

For casts or dressings below the elbow, on older children.

- Seals just above the elbow
- Slim and standard arm options
- Aqua blue junior style
- Suitable for the bath or shower



Age	Weight (kg)	Height	Length (cm)	Circumference (cm) Limb	Cast	Model
8-10 yrs-Slim	Size is typical for age, select as you would for clothing		36	16-22	29	BE810-40
8-10 yrs-Standard			36	22-25	33	BE810-50
11-13 yrs-Slim			43	16-22	29	BE1113-40
11-13 yrs-Standard			43	22-25	33	BE1113-50

Accessories and protection for casts, dressings and PICC Lines

Making life more comfortable



Arm Sleeve for Cast / Dressing / PICC Line

Brighten up your arm cast with a soft lycra cover.

- Helps avoid snagging
- Adult and child sizes
- Elbow / PICC sleeve options - smart and discreet
- Available colours ●●●●●



Adult Leg Sleeve for Cast / Dressing

Keep your leg cast clean with a stylish lycra cover.

- Helps avoid snagging
- Available colours ●



Toe Cozy

Keep toes clean and warm with a cosy fleece cover.

- Adjustable Velcro® straps
- Can be used with dressings too
- Available colours ●●●●●



Oarsome Grips

Extra comfort for crutch handles.

- Helps avoid blisters
- Sold as a pair
- Available colours ●●●●●

Keep clean and dry whilst outside



Adult Foot Cover (Sandal not included)

Keep your cast / dressing clean and dry whilst outside.

- Breathable waterproof fabric
- Adjustable Velcro® straps
- Strong and durable



Adult Mitt / Child Mitt

Keep your cast / dressing clean and dry whilst outside.

- Breathable waterproof fabric
- Adjustable Velcro® straps
- Strong and durable



Adult Wrist Cover

Keep your cast / dressing clean and dry whilst outside with exposed fingers for a better grip.

- Breathable waterproof fabric
- Adjustable Velcro® straps
- Left and right hand options



Waterproof Sock

Rugged outdoor sock for child's lower leg.

- Perfect for the garden or beach
- Waterproof and breathable



Diabetic Sock Collection - Socks that do not constrict blood flow on the lower leg

Prosox

- Stay up without applying pressure
- Welt-free
- Available in cotton or wool
- Anti-slip option available



Gentle Grip Socks

- Non-binding
- Available in black, grey or six vibrant colours as shown
- Extra wide oedema option available




Sock Aid

- Pull on socks with minimum stretching
- Ribbed thumb grips
- Lightweight



Walk-Ease Walking Stick

- Adjustable length
- Handle designed for extra comfort
- Flexible ferrule to maximise contact with the ground
- Three styles



Stick and Stay (Suction Grab Rail)

- Easy to fit
- Ideal for travel
- Suitable for ceramic or glass
- Safety indicator



Helping you whilst you help patients

The bodyblock limb support provides a hygienic multidisciplinary, multipurpose support block for clinical environments.

bodyblock

With four sides creating different heights, the bodyblock is a versatile tool for delivering stable elevation during treatment. Simply flip it for a different height and reduce bending to treat patients. The bodyblock can be used to support the head, legs, heels or arms and can also be used as knee spacer. The bodyblock is lightweight, strong and easy to sterilise.



NO QUIBBLE GUARANTEE

We want users to be totally happy with their LimbO and the service they receive from us. The LimbO is guaranteed to last twelve weeks when used correctly. If the LimbO proves unsuitable for the intended application or fails to last twelve weeks, a refund will be given or replacement will be provided. This does not affect statutory rights.

Please note: The LimbO comes with full instructions and safety information. The LimbO is not intended for use on open or undressed wounds and should not be used in the sea or swimming pools.

Prices correct at October 2022

Here to help you, your team and your patients



Our clinical support team is on hand to help you with:

- Questions about the LimbO
- NHS Supply Chain Order Codes
- Supplying posters, leaflets and literature for your clinic/hospital
- Prescription request forms and FP10 PIP codes
- Samples or demonstrations of the LimbO for your team
- Wholesale/volume purchasing
- Support to assist with CPD
- Details of international re-sellers and distributors

Adult (for quick size selection use height and weight. Use length and circumference points shown on diagram for extra guidance)

Limb Coverage	Full Product Description	Height (cm)	Weight (kg)	Limb Length (cm)	Limb Circumference (cm)	Cast Circumference (cm)	Model	PIP Code	NHS SBS	NHS SC Order Code	Page	Conversion Tables	
Full	Adult Full Arm-Slim	152+	51-70	69	25-29	40 (maximum)	M60L		✓	EAZ 716	4	Stone	Kilos
Full	Adult Full Arm-Medium	152+	70-108	69	30-39	44 (maximum)	M70		✓	EAZ 715	4	8	51
Full	Adult Full Arm-Large	152+	108+	69	39-54	57 (maximum)	M87		✓	EAZ 703	4	8.5	54
Half	Adult Half Arm-Slimmer Adult	152+	Under 51	49	22-25	33 (maximum)	M50A			N/A	4	9	57
Half	Adult Half Arm-Slim to Medium	152+	51-70	53	25-29	40 (maximum)	M60		✓	EAZ 706	4	9.5	60
Half	Adult Half Arm-Medium to Large	152+	70+	53	30-39	44 (maximum)	M67		✓	EAZ 699	4	10	64
Elbow	Adult Elbow-Extra Slim	152+	Under 51	26	22-25	33 (maximum)	M45			N/A	5	10.5	67
Elbow	Adult Elbow-Slim to Medium	152+	51-70	26	25-29	40 (maximum)	M65	386-4303	✓	EAZ 705	5	11	70
Elbow	Adult Elbow-Medium to Large	152+	70-108	26	30-39	44 (maximum)	M75	386-4311	✓	EAZ 711	5	11.5	73
Elbow	Adult Elbow-Extra Large	152+	108+	26	39-54	57 (maximum)	M85			N/A	5	12	76
Full	Adult Hand	Slim to Medium Build		30	22-25	Bandage Only	MITT			N/A	5	12.5	79
Full	Adult Full Leg	165 to 183	64-102	96	52-65	70 (maximum)	M100		✓	EAZ 720	6	13	83
Full	Adult Full Leg-Short	Under 165	64-102	85	52-65	70 (maximum)	M100S			N/A	6	14	89
Full	Adult Full Leg-Long	183+	64-102	103	52-65	70 (maximum)	M100L			N/A	6	14.5	92
Full	Adult Full Leg-Slim	Under 183	Under 64	96	39-54	57 (maximum)	M80LL		✓	EAZ 722	6	15	95
Full	Adult Full Leg-Extra Heavy Leg	152 to 183	102+	96	65-90	94 (maximum)	M200FL		✓	EAZ 721	6	15.5	98
Half	Adult Orthopaedic Boot	For any Orthopaedic Boot please select a full leg Limbo										16	102
Half	Adult Half Leg	165 to 183	64-102	73	39-54	57 (maximum)	M80	320-0045	✓	EAZ 718	6	16.5	105
Half	Adult Half Leg-Short	Under 165	64-102	65	39-54	57 (maximum)	M80S	320-0052	✓	EAZ 707	6	17	108
Half	Adult Half Leg-Long	183+	64-102	83	39-54	57 (maximum)	M80L			N/A	6	17.5	111
Half	Adult Half Leg-Slim	165 to 183	Under 64	73	35-39	44 (maximum)	M76	333-7581	✓	EAZ 710	6	18	114
Half	Adult Half Leg-Slim-Short	Under 165	Under 64	65	35-39	44 (maximum)	M76S	333-7599	✓	EAZ 702	6	18.5	118
Half	Adult Half Leg-Slim-Long	183+	Under 64	83	35-39	44 (maximum)	M76L			N/A	6	19	121
Half	Adult Half Leg-Large	165 to 183	102-140	73	52-57	59 (maximum)	M180	320-0060	✓	EAZ 709	6	19.5	124
Half	Adult Half Leg-Large-Short	Under 165	102-140	65	52-57	59 (maximum)	M180S	320-0078	✓	EAZ 717	6	20	127
Half	Adult Half Leg-Large-Long	183+	102-140	83	52-57	59 (maximum)	M180L			N/A	6	Feet	cm
Half	Adult Half Leg-Extra Heavy Leg	152 to 183	140 +	73	65-90	94 (maximum)	M200HL		✓	EAZ 719	6	4	122
Knee	Adult Knee	165 to 183	64-102	38	39-54	57 (maximum)	M86			N/A	7	4.6	137
Full	Adult Foot-Slim to Medium	Max shoe size 11		20-25	33 (maximum)	M20	416-6476	✓	EAZ 723	7	5	152	
Full	Adult Foot-Medium to Large	Max shoe Size 11		23-34	40 (maximum)	M25	416-6484	✓	EAZ 701	7	5.6	168	
												6	183
												6.6	198
												7	214

Child (for quick size selection use clothing size. Use length and circumference points shown on diagram for extra guidance)

Arm	Full Child Full Arm-2-3 yrs			31	16-22	29 (maximum)	FA23			N/A	9
	Full Child Full Arm-4-5 yrs			43	16-22	29 (maximum)	FA45			N/A	9
	Full Child Full Arm-6-7 yrs			49	22-25	33 (maximum)	FA67			EAZ 712	9
	Full Child Full Arm-8-10 yrs			58	22-25	33 (maximum)	FA810			EAZ 708	9
	Full Child Full Arm-11-13 yrs			58	25-34	40 (maximum)	FA1113			EAZ 713	9
	Half Child Half Arm-8-10 yrs-Slim			36	16-22	29 (maximum)	BE810-40			N/A	9
Elbow	Half Child Half Arm-8-10 yrs-Standard			36	22-25	33 (maximum)	BE810-50			N/A	9
	Half Child Half Arm-11-13 yrs-Slim			43	16-22	29 (maximum)	BE1113-40			N/A	9
	Half Child Half Arm-11-13 yrs-Standard			43	22-25	33 (maximum)	BE1113-50			N/A	9
	Elbow Child Elbow	Under 152	Under 51	12	22-25	33 (maximum)	M45S			N/A	5
Leg	Full Child Full Leg-3 yrs			53	25-34	40 (maximum)	FL3			N/A	8
	Full Child Full Leg-4-5 yrs			58	25-34	40 (maximum)	FL45			N/A	8
	Full Child Full Leg-6-7 yrs			61	35-39	44 (maximum)	FL67			EAZ 714	8
	Full Child Full Leg-8-10 yrs			73	39-54	57 (maximum)	FL810			EAZ 704	8
	Full Child Full Leg-11-13 yrs			83	39-54	57 (maximum)	FL1113			EAZ 700	8
	Half Child Half Leg-8-10 yrs			58	25-34	40 (maximum)	BK810			N/A	8
	Half Child Half Leg-11-13 yrs			61	35-39	44 (maximum)	BK1113			N/A	8

Measuring Guide

